

BD Product Name: **BD Pyxis™ Supply**

Date of Critical or Security Patches: Oct 2021

Abstract: Critical or Security Patches – Oct 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for Oct 2021 will not install on BD Pyxis™ Supply.

The patch includes updates that may adversely affect Pyxis™ Supply. This patch will not be and should not be installed on BD Pyxis™ Supply at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
2021-10 Cumulative Security Update for Internet Explorer 9 for Windows Server 2008 for x86-based systems	This security update resolves vulnerabilities in Internet Explorer. The fixes that are included in this update are also included in the October 2021 Security Monthly Quality Rollup. Installing either this update or the Security Monthly Quality Rollup installs the same fixes.	KB5006671	1
2021-10 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems	This security update includes improvements and fixes that were a part of update KB5005606 (released September 14, 2021) and addresses the following issues: Addresses an issue in which an Internet print server cannot package the driver to send to the client. Implements a Group Policy setting for the following registry value: Registry location: HKEY_LOCAL_MACHINE\Software\Policies\Microsoft\Windows NT\Printers\PointAndPrint Value name: RestrictDriverInstallationToAdministrators	KB5006736	1



	Value data: 1 Adds the ability to configure period or dot (.) delimited IP addresses interchangeably with fully qualified host names in the following Group Policy settings: Package Point and Print - Approved Servers Point and Print Restrictions Addresses an issue in which Security Account Manager (SAM) events are not displayed properly in the Event Viewer.		
2021-10 Security Only Quality Update for Windows Server 2008 for x86-based Systems	This update contains miscellaneous security improvements to internal OS functionality. No additional issues were documented for this release.	KB5006715	1
2021-10 Servicing Stack Update for Windows Server 2008 for x86-based Systems	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. Servicing stack updates (SSU) make sure that you have a robust and reliable servicing stack so that your devices can receive and install Microsoft updates.	KB5006750	1
2021-10 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 and 4.6.2 for Windows Server 2008 SP2	There are no new security improvements in this release. This update is cumulative and contains all previously released security improvements.	KB5006764	1

Note:

1. Due to environment issue patch was unable to be validated this month.

BD Product Name: **BD Pyxis™ CII Safe ES**

Date of Critical or Security Patches: October 2021

Abstract: Critical or Security Patches – October 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2021. The patches include updates that may adversely affect BD Pyxis™ CII Safe ES. These patches will not be and should not be installed on BD Pyxis™ CII Safe ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

Patch Name	Description	Patch ID	Notes
SQL Server 2016 Service Pack 2 Cumulative Update (CU) 16	CU16 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU16 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU16 level.	KB5000645	1
SQL Server 2016 Service Pack 2 Cumulative Update (CU) 15	CU15 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU15 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU15 level. For customers in need of additional installation options, please visit the Microsoft Download Center to download the latest Cumulative Update	KB4577775	1



	(https://support.microsoft.com/en-us/kb/957826)options, please visit the Microsoft Download Center to download the latest Cumulative Update (https://support.microsoft.com/en-us/kb/957826)		
--	---	--	--

Note:

1. After installation of this KB. DB is going into recovery mode. Bug raised for same. Bug Id: 1171472