

Security Patches: BD Pyxis™ Med Station ES

April 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for April 2022 will not install on BD Pyxis™ Med Station ES.

The patch includes updates that may adversely affect Pyxis™ Med Station ES. This patch will not be and should not be installed on BD Pyxis™ Med Station ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

Microsoft® patches

Patch name	Description	Patch ID	Notes
Windows Malicious Software Removal Tool x64 - v5.100	After the download, this tool runs one time to check your computer for infection by specific, prevalent malicious software (including Blaster, Sasser, and Mydoom) and helps remove any infection that is found. If an infection is found, the tool will display a status report the next time that you start your computer. A new version of the tool will be offered every month. If you want to manually run the tool on your computer, you can download a copy from the Microsoft Download Center, or you can run an online version from microsoft.com. This tool is not a replacement for an antivirus product. To help protect your computer, you should use an antivirus product.	KB890830	1
2022-04 Cumulative Update for .NET Framework 4.8 for Windows 10 Version 1607 for x64	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB5011218	1

2022-04 Cumulative Update
for Windows 10 Version
1607 for x64-based
Systems.

This security update includes quality improvements. Key changes include: Addresses a heap leak in PacRequestorEnforcement that degrades the performance of a domain controller. Addresses an issue that affects the Key Distribution Center (KDC) Proxy. The KDC Proxy cannot properly obtain Kerberos tickets for signing in to Key Trust Windows Hello for Business. Addresses an issue that logs Event ID 37 during certain password change scenarios, including failover cluster name object (CNO) or virtual computer object (VCO) password changes. Addresses an issue that causes a Denial of Service vulnerability on Cluster Shared Volumes (CSV). For more information, see CVE-2020-26784. Addresses an issue that prevents you from changing a password that has expired when you sign in to a Windows device.

[KB5012596](#)

1

Note:

1.Medes application is not getting launch on Win10stn14chd12 of Med ES Server 2012 R2-1.4.4.1

Bug Id: 1269876

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

[bd.com](#)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.



Security Patches: BD Pyxis™ Supply

April 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for April 2022 will not install on BD Pyxis™ Supply.

The patch includes updates that may adversely affect Pyxis™ Supply. This patch will not be and should not be installed on BD Pyxis™ Supply at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

Microsoft® patches

Patch name	Description	Patch ID	Notes
2022-04 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems (KB5012658)	<p>This security update includes improvements and fixes that were a part of update KB5011534 (released March 8, 2022) and addresses the following issues:</p> <ul style="list-style-type: none">• Addresses a memory leak that was introduced by the PacRequestorEnforcement registry key in the November 2021 Cumulative Update that causes a decrease in performance on domain controllers.• Addresses an issue in which Event ID 37 might be logged during certain password change scenarios.• Addresses an issue in which domain joins may fail in environments that use disjoint DNS hostnames.• Addresses an issue that prevents you from changing a password that has expired when you sign in to a Windows device.	KB5012658	1

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.





2022-04 Security Only
Quality Update for Windows
Server 2008 for x86-based
Systems (KB5012632)

This security update includes quality improvements.
Key changes include:

[KB5012632](#)

1

- Addresses a memory leak that was introduced by the PacRequestorEnforcement registry key in the November 2021 Cumulative Update that causes a decrease in performance on domain controllers.
- Addresses an issue in which Event ID 37 might be logged during certain password change scenarios.
- Addresses an issue in which domain joins may fail in environments that use disjoint DNS hostnames.

Windows Malicious Software
Removal Tool x64 - v5.100

The Windows Malicious Software Removal Tool (MSRT) helps remove malicious software from computers. Microsoft releases the MSRT on a monthly cadence as part of Windows Update or as a standalone tool. Use this tool to find and remove specific prevalent threats and reverse the changes they have made.

[KB890830](#)

1

Note:

1. Due to environment issue patch was unable to be validated this month.

BD, Franklin Lakes, NJ, 07417, U.S.
201.847.6800

bd.com

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.