

BD Product Name: **BD Pyxis™ CII Safe ES**

Date of Critical or Security Patches: May 2021

Abstract: Critical or Security Patches – May 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for May 2021. The patches include updates that may adversely affect BD Pyxis™ CII Safe ES. These patches will not be and should not be installed on BD Pyxis™ CII Safe ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

| Patch Name | Description | Patch ID | Notes |
|--|---|---------------------------|-------|
| Security Update for SQL Server 2016 Service Pack 2 CU | This security update fixes the following issue: KB4583468 - Microsoft SQL Server elevation of privilege vulnerability | KB4583461 | 1 |
| SQL Server 2016 Service Pack 2 Cumulative Update (CU) 16 | CU16 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU16 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU16 level. | KB5000645 | 1 |
| SQL Server 2016 Service Pack 2 Cumulative Update (CU) 15 | CU15 for SQL Server 2016 Service Pack 2 upgraded all SQL Server 2016 Service Pack 2 instances and components installed through the SQL Server setup. CU15 can upgrade all editions and servicing levels of SQL Server 2016 Service Pack 2 to the CU15 level. For customers in need of additional installation options, please visit the Microsoft Download Center | KB4577775 | 1 |



| | | | |
|--|--|--|--|
| | to download the latest Cumulative Update (https://support.microsoft.com/en-us/kb/957826)options, please visit the Microsoft Download Center to download the latest Cumulative Update (https://support.microsoft.com/en-us/kb/957826) | | |
|--|--|--|--|

Note:

1. After installation of this KB. DB is going into recovery mode. Bug raised for same. Bug Id: 1171472

BD Product Name: **BD Pyxis™ MedStation™ ES**

Date of Critical or Security Patches: May 2021

Abstract: Critical or Security Patches – May 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for May 2021. The patches include updates that may adversely affect BD Pyxis™ MedStation™ ES. These patches will not be and should not be installed on BD Pyxis™ MedStation™ ES at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering.

| Patch Name | Description | Patch ID | Notes |
|---|--|---------------------------|-------|
| 2021-05 Cumulative Update for Windows 10 Version 1809 for x64-based Systems | A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system. | KB5003171 | None |
| 2021-04 Cumulative Update for Windows 10 Version 1809 for x64-based Systems | A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system. | KB5001342 | None |



BD Product Name: **BD Pyxis™ Supply**

Date of Critical or Security Patches: May 2021

Abstract: Critical or Security Patches – May 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft® that have been identified as critical or security related for May 2021 will not install on BD Pyxis™ Supply.

The patch includes updates that may adversely affect Pyxis™ Supply. This patch will not be and should not be installed on BD Pyxis™ Supply at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only apply to the latest supported version(s) of the BD offering.

| Patch Name | Description | Patch ID | Notes |
|--|--|---------------------------|-------|
| 2021-05 Cumulative Security Update for Internet Explorer 9 for Windows Server 2008 for x86-based systems | This security update resolves vulnerabilities in Internet Explorer. | KB5003165 | 1 |
| 2021-05 Security Only Quality Update for Windows Server 2008 for x86-based Systems | This security update includes quality improvements. Key changes include: Security updates to Windows App Platform and Framework. | KB5003225 | 1 |
| 2021-05 Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems | This security update includes improvements and fixes that were a part of update KB5001389 (released April 13, 2021) and addresses the following issues: Security updates to Windows App Platform and Framework, and Microsoft Scripting Engine. | KB5003210 | 1 |
| 2021-05 Security and Quality Rollup for .NET Framework 2.0, 3.0, 4.5.2, 4.6 for Windows Server 2008 SP2 | The update addresses the vulnerability by correcting how the .NET Framework handles objects in memory. | KB5001882 | 1 |

Note:

1. Due to environment issue patch was unable to be validated this month.

