

# Security Patches: BD Pyxis™ CIISafe

Oct 2022

BD has identified patches from Microsoft® that have been identified as critical or security related for October 2022 will not install on BD Pyxis™ CIISafe. The patch includes updates that may adversely affect Pyxis™ CIISafe. This patch will not be and should not be installed on BD Pyxis™ CIISafe at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: This patch only applies to the latest supported version(s) of the BD offering.

## Microsoft® patches

Patch name	Description	Patch ID	Notes
2022-09 Cumulative Update for Windows 10 Version 1607 for x64-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	<a href="#">KB5017305</a>	1
2022-09 Servicing Stack Update for Windows 10 Version 1607 for x64-based Systems	Install this update to resolve issues in Windows. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article for more information. After you install this item, you may have to restart your computer.	<a href="#">KB5017396</a>	1



2022-10 Cumulative Update  
for Windows 10 Version  
1607 for x64-based Systems

A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.

[KB5018411](#)

1

Windows Malicious Software  
Removal Tool x64 - v5.106

After the download, this tool runs one time to check your computer for infection by specific, prevalent malicious software (including Blaster, Sasser, and Mydoom) and helps remove any infection that is found. If an infection is found, the tool will display a status report the next time that you start your computer. A new version of the tool will be offered every month. If you want to manually run the tool on your computer, you can download a copy from the Microsoft Download Center, or you can run an online version from [microsoft.com](https://microsoft.com). This tool is not a replacement for an antivirus product. To help protect your computer, you should use an antivirus product.

[KB890830](#)

1

Note:

1. Patches not getting downloaded and installed on Win 10 Physical CIISafe Legacy devcie.

BD, Franklin Lakes, NJ, 07417, U.S.  
201.847.6800

[bd.com](https://www.bd.com)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2022 BD. All rights reserved.

