

BD Pyxis™ MedStation™ 4000

December 2019

Critical and/or Security Patches – December 2019

BD has identified patches from Microsoft® that have been identified as critical or security related for December 2019.

The patches include updates that may adversely affect MedStation™4000. These patches will not be and should not be installed on MedStation™4000 at this time.

Customers that maintain patches independent of BD automated delivery should ensure these actions are not performed as the acting responsible entity in order to maintain the correct security posture of the system(s).

Note: These patches only apply to the latest supported version(s) of the BD offering

Patch Name	Patch Description	Patch ID	Notes
Security Monthly Quality Rollup for Windows Server 2008 for x86-based Systems	A security issue has been identified in a Microsoft software product that could affect your system. You can help protect your system by installing this update from Microsoft. For a complete listing of the issues that are included in this update, see the associated Microsoft Knowledge Base article. After you install this update, you may have to restart your system.	KB4530695	1

