

## BD Synapsys™

### Date of Critical or Security Patches: December 2021

Critical Security Patches: September 2021 - October 2021

### Microsoft® & Third-Party Patches for Windows Server 2016

BD has identified patches from Microsoft® that have been identified as critical or important security related updates thru October 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Patch Description	Patch ID	Notes
Cumulative Security Update	Updates security for the Windows operating system.	<a href="#">KB5006669</a>	
Servicing Stack Update	This update makes quality improvements to the servicing stack, which is the component that installs Windows updates.	<a href="#">KB5005698</a>	1

#### Notes:

1. Servicing Stack update KB5005698 must be installed prior to other updates.

