

BD Synapsys™

Date of Critical or Security Patches: September 2021

Critical Security Patches: January 2021-July 19, 2021

Microsoft® & Third-Party Patches for Windows Server 2016

BD has identified patches from Microsoft® that have been identified as critical or important security related updates thru March 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

| Patch Name | Patch Description | Patch ID | Notes |
|-------------------------------------|--|---------------------------|-------|
| Security Update for Secure Boot DBX | This security update makes improvements to Secure Boot DBX | KB4535680 | |
| Cumulative Security Update | <ul style="list-style-type: none">▪ Updates to remove Adobe Flash from your device.▪ Updates for verifying usernames and passwords.▪ Updates to improve security when Windows performs basic operations. | KB5004238 | |
| SQL Server 2016 Update | This update contains fixes that were released after the initial release of SQL Server 2016 SP2. | KB5001092 | |
| Servicing Stack Update | This update makes quality improvements to the servicing stack, which is the component that installs Windows updates. | KB5001402 | 1 |

Notes:

Servicing Stack update KB5001078 must be installed prior to other updates.

