

BD Product Name: **Identity provider management (IDM)**

Date of Critical or Security Patches: December 2021

Abstract: Critical or Security Patches – December 2021

Microsoft® & Third-Party Patches

BD has identified patches from Microsoft that have been identified as critical or security related for December 2021. These patches were not found to adversely affect BD products and will be applied according to customers' service agreement.

Customers that maintain patches independent of BD automated delivery should ensure the validated patches are installed on their BD systems as the acting responsible entity in order to maintain the correct security posture of the system(s).

Patch Name	Description	Patch ID	Notes
2021-12 Security Only Quality Update for Windows Server 2012 R2 for x64-based System	<p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none">▪ Update to support the cancellation of daylight savings time (DST) for 2021 for the Republic of Fiji.▪ Addresses a known issue that causes error codes 0x000006e4, 0x0000007c, or 0x00000709 when connecting to a remote printer that is shared on a Windows print server.▪ Addresses a known issue that might prevent apps, such as Kaspersky apps, from opening after you attempt to repair or update the apps using the Microsoft Installer (MSI). <p>For more information about the resolved security vulnerabilities, please refer to the new Security Update</p>	KB5008285	N/A

	Guide website and the December 2021 Security Updates.		
2021-12 Security Only Quality Update for Windows Server 2012 R2 for x64-based System	<p>This security update includes improvements and fixes that were a part of update KB5007247 (released November 9, 2021) and addresses the following issue:</p> <ul style="list-style-type: none"> ▪ Update to support the cancellation of daylight savings time (DST) for 2021 for the Republic of Fiji. ▪ Addresses a known issue that causes error codes 0x000006e4, 0x0000007c, or 0x00000709 when connecting to a remote printer that is shared on a Windows print server. ▪ Addresses a known issue that might prevent apps, such as Kaspersky apps, from opening after you attempt to repair or update the apps using the Microsoft Installer (MSI). 	KB5008263	N/A
Windows Malicious Software Removal Tool x64 - v5.96	<p>The Windows Malicious Software Removal Tool (MSRT) helps remove malicious software from computers that are running any of the following operating systems:</p> <ul style="list-style-type: none"> ▪ Windows 10 ▪ Windows Server 2019 ▪ Windows Server 2016 ▪ Windows 8.1 	KB890830	N/A

	<ul style="list-style-type: none"> ▪ Windows Server 2012 R2 ▪ Windows Server 2012 R2 ▪ Windows Server 2008 R2 ▪ Windows 7 ▪ Windows Server 2008 		
2021-12 Cumulative Update for Windows Server 2016 for x64-based Systems	<p>This security update includes quality improvements. Key changes include:</p> <ul style="list-style-type: none"> ▪ Adds support for the cancellation of daylight savings time for the Republic of Fiji for 2021. ▪ Addresses a known issue that causes error codes 0x000006e4, 0x0000007c, or 0x00000709 when connecting to a remote printer that is shared on a Windows print server. ▪ Addresses a known issue that might prevent apps, such as Kaspersky apps, from opening after you attempt to repair or update the apps using the Microsoft Installer (MSI). 	KB5008207	N/A