




Arctic Gel Pads Par Level Order Sheet:

Hospital name:


Arctic Gel™ Kit Sizes	Order Number	Internal Stock Number	Par Level	Current Stock	Need to Order
Large	#317-09-02				
Medium	#317-07-02				
Small	#317-05-02				
X-Small	#317-03-02				
XX Small	#317-02-02				
Universals	#317-00-04				
Ped Universal	#318-01-04				
Neonatal Pads	#318-02-02				

- ArcticGel™ Pads come in “Kits”. Each Kit includes 4 pads total, 2 chest pads and 2 legs pads to cover one patient
- Kits are packaged in a Case. Each Case has two “Kits”
- Universal pads are packaged individually and separately as needed for additional coverage or to replace leg pads as needed. There are 4 Universal pads per case.
- Neonatal Pad is a single pad and comes 2 neonatal pads per case
- Expiration date is 24 months







Neonatal
1.8-4.5kg
(4.0-9.9 lbs)
318-02-02




1 Small Universal
2.5-5kg
(5.5-10 lbs)
318-01-04




2 Small Universal
5-10kg
(10-22 lbs)
318-01-04




3 Small Universal
10-16kg
(22-35 lbs)
318-01-04




XX Small
16-30kg
(35-66lbs)
317-02-02




X Small
31-45kg
(67-98 lbs)
317-03-02




Small
44-57kg
(98-125 lbs)
317-05-02



Medium
57-73kg
(126-160 lbs)
317-07-02



Large
73-100kg
(161-220 lbs)
317-09-02



Large + Universal
100+kg
(220+ lbs)
317-00-04 (Universal)

Follow-up notes:

TM Name:

CM Name:

TM Email:

CM Email:

TM Phone #:

CM Phone #:

Please consult product inserts and labels for any indications, contraindications, hazards, warnings, cautions and directions for use. Federal Law (USA) restricts this device to sale by or on the order of a physician.

Medivance, Inc.
321 South Taylor Ave., Suite 200
Louisville, Colorado 80027 USA

bd.com

Customer Service Phone: 1.844.823.5433
Customer Service Fax: 1.800.852.1339
Email: arcticsun@bd.com

BD, the BD logo, Arctic Sun, and ArcticGel are trademarks of Becton, Dickinson and Company and its affiliates.
©2022 BD. All Rights Reserved. 2211-01 BD-77237

