


# BD PhaSeal™ Optima Closed System Drug Transfer Device (CSTD)

## Connection technique for administration

Airtight and leakproof, and built on the legacy of the BD PhaSeal™ CSTD, the BD PhaSeal™ Optima system mechanically prohibits the transfer of environmental contaminants into the system and the escape of drug vapor concentrations outside the system, thereby minimizing individual and environmental exposure to drug vapor, aerosols and spills.

### Engage



### Administer



### Disengage



- Make sure you hear the audible click when engaging the Connector/Injector
- Remember to check your connections