

Texium® IV administration sets

Close your system using pre-assembled Texium® administration sets

Texium IV administrations sets provide a pre-assembled three-in-one solution, combining *three products in one* to help optimize safety, simplify workflow and minimize costs associated with reducing exposure to hazardous drugs.



Three-in-one solution

Bonded Texium administration sets provide three products in one:

- Alaris® SmartSite® dedicated primary pump or gravity tubing available in low-sorbing or conventional DEHP-Free tubing
- SmartSite drug transfer port fully integrated with the drip chamber
- Texium closed male luer pre-bonded to the end of intravenous (IV) tubing

Additionally, bonded Texium administration sets come with a priming cap to help assist with the priming of the set.

Benefits

- No component assembly, streamlining workflow
- Permanently bonded, reducing the risk of disconnect, leakage and exposure
- Needle-free luer locking, complying with the National Institute for Occupational Safety and Health (NIOSH)/Occupational Safety and Health Administration (OSHA) standards and improving nursing and patient safety
- Fewer parts, decreasing cost and hazardous waste
- Requires no change to nursing protocol and workflow
- Passive locking mechanism, helps reduce leakage and exposure and preventing accidental leaks
- Seamless connection with SmartSite y-site valve integrated with Alaris pump tubing
- Available with low-sorbing tubing



CareFusion
San Diego, CA

carefusion.com

