



BD Pyxis™ Full-Height CUBIE™ Pocket

What Fits

Medication type	Packaging*	Pocket size**	Quantity
Nystatin	5 mL unit dose (<i>small</i>)	2x1	18–20
Nystatin	5 mL unit dose (<i>large</i>)	(105 x 76 x 55)	8
Albuterol sulfate inhalation aerosol	90 mcg per inhalation		4 (<i>in pouch packaging</i>)
Hydromorphone	5 mL syringe		12–15
Nystatin	5 mL unit dose (<i>small</i>)	2x2	40+
Propofol	100 mL vials (<i>bottles fit lying down</i>)	(105 x 69 x 12)	2
Propofol	20 mL vials		12
Pharmacy-compounded product	20 mL syringe		8–9
Ibuprofen suspension	4 oz bottle		4
Tylenol with codeine	15 mL unit dose		24–25
Fentanyl/Morphine drip	100 mL (<i>in 150 mL bag</i>)		4–5
Heparin	500 mL bag	2x3	1
Lovenox	Syringe	(105 x 69 x 194)	10–12
Magnesium sulfate	100 mL bag		8
Propofol	20 mL vials		18–19
Tylenol with codeine	30 mL unit dose		14
Patient-controlled analgesia (PCA)	30 mL syringe		10–12
Heparin	500 mL bag		2
PCA	50–60 mL syringe		5–6
Propofol	100 mL vials (<i>bottles fit lying down</i>)	2x5	6
Tylenol with codeine	30 mL unit dose	(105 x 69 x 337)	23
Vancomycin vial pre-attached to IV bag	250 mL bag + vial		3

*Medication packaging may vary by manufacturer. Represented quantity is an approximation.

**Minimum internal H x D x W dimensions in mm.



BD Pyxis™ Half-Height CUBIE™ Pocket

What Fits

Medication type	Pocket size**		
	1x1 (38 x 86 x 53)	1x2 (38 x 86 x 117)	1x3 (38 x 86 x 183)
Small tablet—manufacturer UD	30/tear into 6s/4s	80s/10s	Combine***
Large tablet—manufacturer UD	20/tear into 6s/4s	50s/10s	Combine
Small capsule—manufacturer UD	24s/6s and some 2s/1s	80s/10s	Combine
Large capsule—manufacturer UD	20/tear into 6s/4s	50s/10s	Combine
Small tablet/cap—small blister pack	20 singles	40 singles	60 singles
Large tablet/cap—large blister pack	15 singles	30 singles	45 singles
Patches	Small patches < 15	20+	Combine
Vials 2 mL	15 singles	30 singles	45 singles
Vials 4 mL	8 singles	16 singles	24 singles
Vials 10 mL	3 singles	6 singles	9 singles
Vials 20 mL	2 singles	4 singles	6 singles
Vials 30 mL	1 single	3 singles	4 singles
Mix-o-vials (e.g., Solu-Medrol®)	3 singles	9 singles	12 singles
Ampoules 1 mL	20 singles	40 singles	60 singles
Ampoules 2 mL	16–18 singles	36–40 singles	Combine
Ampoules 5 mL	8 singles	16 singles	24 singles
Syringes	0	0	Approx 10 pre-filled
Advantage vials	1 single	4 singles	5 singles
I-Pak	0	30 singles	Approx 50 singles
Vials 1 g antibiotics	1, 2 or 3, depending on size	4–7	Combine
Vials 2 g antibiotics	1–2	3–4	Combine
6 dram prescription vials with flip tops only	1	1–2	1–3
10 dram prescription vials with flip tops only	1	1	1–2

*Medication packaging may vary by manufacturer. Represented quantity is an approximation. Capacity information above applies to all styles of BD Pyxis CUBIE hardware.

**Minimum internal H x D x W dimensions in mms.

***Combine represents the combined quantity that fits in a 1x1 and 1x2 pocket size.

bd.com

