

# When medical emergencies arise, every minute counts

## BD Vacutainer® Rapid Serum Tube with Hemogard™ Closure



### Improved turnaround time

- Clots in only 5 minutes!
- 3 minute spin\*
- May reduce laboratory turnaround time by up to 32 minutes\*\*

### Standardize chemistry stat testing to 1 tube

- The Speed of Plasma, the Quality of Serum
  - 1 draw, 1 tube

### Improve specimen quality

- Improving patient care and ED productivity
  - 59% reduction in hemolysis\*\*\*
  - 59% reduction in fibrin strands\*\*\*
  - Significant reduction in Red Blood Cell contamination\*\*\*



## Improve turnaround time; and potential reduction in redraws

- Improved patient care, ED throughput
- ↳ Improved patient satisfaction
- ↳ Increased ED revenue

\* At 4000 x g (RCF) for 3 minutes

\*\* Compared to current serum gel tubes

\*\*\* Compared to control serum gel tubes, based on a meta-analysis of lots from multiple studies.

Data on File



Helping all people  
live healthy lives


# BD Vacutainer® Rapid Serum Tube

with Hemogard™ Closure

Feature	Benefit
① BD Hemogard™ Closure	<b>Improves Health Care Worker (HCW) safety by reducing:</b> <ul style="list-style-type: none"> <li>• Contact with blood on the stopper</li> <li>• Blood splatter during de-capping</li> </ul> Compatible with major laboratory automation systems
② V-Notch™ Label	• Helps align secondary label
③ Thrombin based clot activator	• Promotes rapid clotting
④ Separator	<ul style="list-style-type: none"> <li>• DEHP-free gel</li> <li>• Extended analyte stability</li> <li>• Instrument automation compatible</li> </ul>



Ref #: 368774



BD Vacutainer® Rapid Serum Tube is **ready for analysis in only 8 minutes.**  
**32 valuable minutes saved** with RST.

**For more information, call or visit us online:**

BD Global Technical Services at **1.800.631.0174**  
 or e-mail [www.bd.com/vacutainer/contact](http://www.bd.com/vacutainer/contact)  
 BD Customer Service at **1.888.237.2762**  
 or visit us anytime online at [www.bd.com/vacutainer](http://www.bd.com/vacutainer)



**BD Diagnostics**  
 Preanalytical Systems  
 1 Becton Drive  
 Franklin Lakes, NJ 07417  
[www.bd.com/vacutainer](http://www.bd.com/vacutainer)