

When every minute counts



BD Vacutainer® Rapid Serum Tube with Hemogard™ Closure

Improved turnaround time

- Clots in only 5 minutes!
- 3 minute spin*
- May reduce laboratory turnaround time by up to 32 minutes**

Improved lab efficiency

- Fewer tubes drawn, fewer reworks
 - Enhanced streamlining of workflow

Improve specimen quality

- Improved lab productivity and patient care
 - 59% reduction in hemolysis***
 - 59% reduction in fibrin strands***
 - Significant reduction in Red Blood Cell contamination***

Standardize chemistry stat testing to 1 tube

- The Speed of Plasma, the Quality of Serum
 - 1 draw, 1 tube



Faster turnaround time → Impact decision making
→ Improved patient care → Improved patient satisfaction

* At 4000 x g (RCF) for 3 minutes

** Compared to current serum gel tubes

*** Compared to control serum gel tubes, based on a meta-analysis of lots from multiple studies.

Data on File



Helping all people
live healthy lives


BD Vacutainer® Rapid Serum Tube

with Hemogard™ Closure

Feature	Benefit
1 BD Hemogard™ Closure	Improves Health Care Worker (HCW) safety by reducing: <ul style="list-style-type: none"> • Contact with blood on the stopper • Blood splatter during de-capping Compatible with major laboratory automation systems
2 V-Notch™ Label	• Helps align secondary label
3 Thrombin based clot activator	• Promotes rapid clotting
4 Separator	<ul style="list-style-type: none"> • DEHP-free gel • Extended analyte stability • Instrument automation compatible



Ref #: 368774



BD Vacutainer® Rapid Serum Tube is **ready for analysis in only 8 minutes.**
32 valuable minutes saved with RST.

For more information, call or visit us online:

BD Global Technical Services at **1.800.631.0174**
 or e-mail www.bd.com/vacutainer/contact
 BD Customer Service at **1.888.237.2762**
 or visit us anytime online at www.bd.com/vacutainer



BD Diagnostics
 Preanalytical Systems
 1 Becton Drive
 Franklin Lakes, NJ 07417
www.bd.com/vacutainer