



Because your patients trust
you to be your best

BD Anesthesia Solutions



About BD

Anesthesia Solutions

A tradition of excellence in regional anesthesia

BD has produced regional anesthesia products for most of our 118-year history. We have the experience, expertise and technology necessary to help you in your practice.

A full range of innovative choices

BD offers a wide range of regional anesthesia devices and components to meet your needs. We have a proven track record of innovation in needle technology and a market-driven product development process.

Cost-efficient programs backed by specialists

- The breadth and quality of BD brand components provide both clinical and economic value
- The BD Regional Anesthesia Program provides a wide range of Regional Anesthesia Trays
- BD educational materials help provide a foundation for procedural training
- BD sales consultants can help provide both technical and clinical solutions

Dealing with the issues that affect you most

Read on to learn about the products and programs BD offers. Call your BD sales consultant for more information or to set a clinical assessment



BD™ Whitacre Pencil Point Spinal Needles



BD Whitacre Pencil Point

Available in gauges from 22 to 27 and in lengths from 3 1/2" to 5"

BD Whitacre Spinal needles

Standard length

Description	Order no.
22 G x 3 1/2"	405010
24 G x 3 1/2" high flow	405133
25 G x 3 1/2" high flow	405138
27 G x 3 1/2" high flow	405079

Long length

Description	Order no.
25 G x 4 11/16"	409442
25 G x 5" high flow*	405140
27 G x 4 11/16" high flow	409443
27 G x 5" high flow*	405144

Spinal introducers

Description	Order no.
20 G x 1 1/4" (use with 25 G and 27 G needle or smaller)	405160

Available in BD Regional Anesthesia Trays

*Not available in Trays

BD Quincke Spinal Needles

BD Quincke Spinal needles

Standard length	
Description	Order no.
18 G x 3" *	405174
18 G x 3 1/2"	405184
20 G x 2 1/2"	405071
20 G x 3" *	405172
20 G x 3 1/2"	405182
22 G x 2 1/2" *	405074
22 G x 3"	405171
22 G x 3 1/2"	405181
23 G x 3 1/2"	400106
25 G x 3" *	405170
25 G x 3 1/2"	405180
26 G x 3 1/2"	405164
27 G x 3 1/2"	405081

Long length	
Description	Order no.
18 G x 6" *	408360
20 G x 6" *	405211
22 G x 5" *	405148
22 G x 7" *	405149
25 G x 4 11/16" *	405234

Short length	
Description	Order no.
20 G x 1 1/2" *	405162
22 G x 1 1/2" *	405161
25 G x 1" *	405073
25 G x 2" *	405078

Spinal introducers	
Description	Order no.
20 G x 1 1/4" (use with 25 G and 27 G needle or smaller)	405160

Available in BD Regional Anesthesia Trays

*Not available in Trays



BD Quincke Point

Available in gauges from 18 to 27 and in lengths from 1" to 7"

BD™ Perisafe Epidural Needles



BD Tactile Tuohy Point



BD Modified Tuohy Point

BD Perisafe™ Tactile Epidural needles

The BD™ Perisafe Tactile Epidural needle is specifically designed for those clinicians who prefer an epidural needle with tactile feedback.

The BD™ Tactile Tuohy Epidural needle tip has a primary and secondary bevel. The secondary bevel design alters the geometry of the needle and provides tactile feedback while maintaining the sharpness to penetrate underlying tissues.

BD Modified Tuohy Point needles

The BD Perisafe™ Modified Tuohy Point Epidural needle with point geometry employs a primary bevel. It is designed for those clinicians with a focus on penetration characteristics.

BD Perisafe™ Epidural needles

Weiss (fixed wing) needles

Description	Order no.
17 G x 3 1/2" tactile	405528
18 G x 3 1/2" tactile	405529
17 G x 3 1/2" modified	405191
17 G x 5" modified	405193
18 G x 3 1/2" modified	405188
18 G x 5" modified	405190
20 G x 2" modified	405185

Tuohy (non-winged) needles

17 G x 3 1/2" tactile	405530
18 G x 3 1/2" tactile	405531
17 G x 3 1/2" modified	405092
18 G x 3 1/2" modified	405018
20 G x 3 1/2" modified	405028

Available in BD Regional Anesthesia Trays

BD Epidural Products

BD Perisafe™ catheter with catheter connector

BD Perisafe™ catheter and connector is available in a sterile pouch with a closed-end multi-port tip configuration, thread assists and connector.



BD Perisafe™ nylon epidural catheters with connectors


Description	Order no.
20 G x 36" closed end for use with 18 G x 3 1/2" needle	400507

Available in [BD Regional Anesthesia Trays](#)

Spring wound epidural catheter

- Choice of soft, rounded, closed tip with multiple ports, or open tip available in a variety of tray configurations
- Large, 1" proximal and 2" open series of coils at the distal tip to allow visualization of fluids
- Blended nylon over stainless steel wound spring construction with circumferential centimeter markings

Only Available in [BD Regional Anesthesia Trays](#)



BD Epilor™ loss of resistance syringe

- A double ridge syringe stopper assembly that provides a secure seal to the syringe barrel, allowing both aspiration and injection of air and / or saline
- A positive plunger rod stop to prevent accidental disassembly

BD Epilor™ loss of resistance syringe


Description	Order no.
7 mL Luer-Lok™, plastic	405291
7 mL luer-slip, plastic	405292

Available in [BD Regional Anesthesia Trays](#)



BD Epilor Plastic Syringe


BD Combined Spinal / Epidural System




12 mm extension when hubs are fully engaged



BD Durasafe Plus™ Epidural Lock CSE Needle Set



Lock release button
Variable extension adapter



Unlock Matched luer-slip hubs

Adjusts spinal needle length from epidural needle

BD Durasafe™ CSE Needle Set

- CSE matched set system with BD™ Whitacre high flow pencil point needle
- BD Durasafe™ matched sets designed to allow up to 12 mm extension

BD Durasafe™ Plus Epidural Lock CSE Needle Set

BD Durasafe™ needle sets	
Description	Order no.
25 G x 4 11/16" Whitacre / 17 G x 3 1/2" Weiss needle	405139
27 G x 4 11/16" Whitacre / 18 G x 3 1/2" Weiss needle	405452
27 G x 4 11/16" Whitacre / 17 G x 3 1/2" Weiss needle	405453

- Epidural lock needle permits the spinal needle to lock at any extension
- Epidural lock needle helps ensure securement of the spinal needle after dural puncture

BD Durasafe™ Plus Epidural Lock CSE Needle Set

BD Durasafe Plus™ Epidural Lock CSE Needle Set	
Description	Order no.
27 G BD Whitacre, 17 G Weiss Locking needle, Closed end 19 G catheter, BD Epilor™ plastic LOR	400286
27 G BD Whitacre, 18 G Weiss Locking needle, Closed end 20 G catheter, BD Epilor™ plastic LOR	400287

Available in BD Regional Anesthesia Trays

BD Peripheral Block Needles

Nerve block needles

- Available in insulated and uninsulated configuration for surface and deep tissue blocks
- Block bevel enables skin penetration with short geometry for medication delivery
- Insulated BD Stimex™ needles with block bevel available in 2 1/8" length

Nerve block needles	
Description	Order no.
22 G x 2 1/8" (insulated, block bevel)	404204
22 G x 1 1/2" (uninsulated, block bevel)	408348

BD™ Ultra-Vue needles

- Echogenic general use procedure needles
- Designed for use with ultrasound guidance

BD™ Echogenic procedure needles (for use with ultrasound)	
Description	Order no.
BD™ Ultra-Vue 22 G x 3 1/2" spinal type point	400727
BD™ Ultra-Vue 20 G x 3 1/2" spinal type point	400730

Specialty needles

- BD Chiba™ needles available with short bevel point geometry offered in 6" and 8" lengths
- Short bevel and hypodermic bevel needles also available in a full range of sizes and gauges

BD Chiba bevel needles	
Description	Order no.
22 G x 6"	408268
22 G x 8"	408269

Short bevel needles	
Description	Order no.
22 G x 1 1/2"	305159
20 G x 1 1/2"	305179
18 G x 1 1/2"	305199




Block bevel



BD™ Ultra-Vue echogenic needle



BD Chiba bevel



Short bevel needle

BD Clinical Specialty Needles

BD™ Breast Lesion Localization needles

Description	Order no.
BD PrecisionGlide™ 20 G x 3 1/2" (8.9 cm) with hook wire	409407
BD PrecisionGlide™ 20 G x 2 1/2" (6.4 cm) with hook wire	409408
BD PrecisionGlide™ 20 G x 1 1/2" (3.8 cm) with hook wire	409409

BD™ Liver, Kidney and Soft Tissue (reusable needles)


Description	Order no.
14 G x 3 3/8" Silverman, Modified Franklin, thin wall	511548
14 G x 4 1/2" Silverman, Modified Franklin, thin wall	511547
16 G x 4" Klatskin, thin wall	511403

BD™ Pleural Biopsy needle (reusable needles)

Description	Order no.
11 G x 3" Cope	511421

BD™ Bone Marrow, Aspiration and Transfusion* (reusable needles)

18 G x 1" Rosenthal	511341
16 G x 1 5/16" Rosenthal	511346
18 G x 1" Osgood	511348
16 G x 1 5/16" Osgood	511349
10 G x 3" Large Silverman	511420
10 G x 4 1/4" Westerman-Jensen	511573



BD Chiba Fine

BD™ Angiography needles (for Seldinger Technique)

Description	Order no.
POTTS-COURNAND needle, thin wall with obturator and cannula extender	
18 G x 2 3/64" (.032" - .038") guide wire	408270
Arterial Needle Outer Blunt Cannula	
18 G x 2 7/8" (.032" - .038") guide wire	408273
Single wall puncture (M / P) procedure needle	
18 G x 2 3/4" one part needle with wings (.032" - .038") guide wire	408295
18 G x 2 3/4" one part needle (.032" - .038") guide wire	408300

BD™ Westcott Fine Needle Aspiration Biopsy

22 G x 3 1/2"	408265
22 G x 5 1/2"	408266
22 G x 7"	408267
20 G x 3 1/2"	408262
20 G x 5 1/2"	408263
20 G x 7"	408264

BD™ Chiba Fine Needle Aspiration Biopsy / Transhepatic Cholangiography needles

Description	Order no.
22 G x 6" Chiba type, thin wall	408268
22 G x 8" Chiba type, thin wall	408269

BD Regional Anesthesia Trays

High quality product, dependable service and a broad portfolio

Value-added and customized solutions

BD is committed to providing solutions for your regional anesthesia practice needs through a broad array of procedural trays. BD trays contain the same recognized and trusted components for spinal, epidural, combined spinal and epidural procedures.

Using our proprietary customer focused needs analysis tools, BD sales consultants can help you:

- Review your product usage and procedural preferences
- Identify needs in components selection and procedure configuration



Page is blank so that next
page is on the back

For a clinical assessment, contact
your local BD sales consultant

BD, Franklin Lakes, NJ, 07417, US

bd.com

