

35	TC10008126	Board, Display ETCO2.
55	TC10008564	Board, Logic
39	TC10008127	Board, Power Supply
20	49000215	Case, Front (Front Case/Keypad Assembly)
85	10013667	Case, Rear (Assembly Nameplate & Labels)
24	321246	Cotter Pin
21	TC10003415	Door Assembly
23	148228-000	Door Assembly, Spring Retainer
22	148130-002	Door Assembly, Torsion Spring
47	148249-100	EtCO2, Cable Assembly
45	148131-001	Frame Mech
40	TC10008320	IUI Connector Bracket
60	146992-001	IUI Connector Bracket Strap
130	147078-100	IUI Connector Kit Left (IUI Connector Left and Gasket)
120	147077-100	IUI Connector Kit Right (IUI Connector Right and Gasket)
90	147080-100	Latch Kit (Consists of: Latch, Latch Support, Leaf Spring and Compression Spring)
32	144424-100	Keypad Backlight LED Assembly 8220/8210
635	10010621	Label, Door
630	148372-000	Label, Microstream
620	148405-000	Label, Patent
16	148140-000	Luer Retainer
844	300928	Nut, Hex (Screw, Luer Retainer)
37	11754144	Oridion Module Kit
17	321245	O-Ring .364ID x .070W 60 DURO
46	TC10003413	PCB Holder (Retainer)
75	142794-000	Rubber, Feet
825	300348	Screw, (Frame)
841	300923	Screw, (IUI Bracket Ground Strap) & (Display Board)
835	320855	Screw, (IUI Connector) & (Rear Case)
845	300345	Screw, (Latch Screws)
842	303822	Screw, (Luer Retainer)
846	305121	Screw, (Oridion Module, PCB Retainer)
843	300347	Screw, (Oridion Module, Stand-Off)
830	806112	Snap Rivet
36	321270	Standoff F/F 3/16 HEX x 3/8 4-40 AL
25	146978-001	Status Indicator Lens
N/A	146980-000	Status Indicator Lens Gasket
N/A	10015532	Technical Service Manual
80	147093-100	Tubing MED Silicon, 10 feet
840	320908	Washer, (Rear Case Screw) & (Oridion Module, Stand-Off Screw Washer)