

3	49000227	Actuator Knob Assembly Kit (Complete Assembly)
255	TC10003065	Actuator Knob Assembly (Component Only)
275	147836-002	Actuator Knob, (Cam Pull Rod)
271	321041	Actuator Knob, (Compression Spring)
265	TC10008324	Actuator Knob, (Flag Shaft Interrupter)
230	147799-001	Actuator Knob, (Plate Slide)
397	320999	Assembly, Cable Dual End Connector
325	TC10006400	Barrel Clamp (Use TC10005065 until inventory depleted)
342	147987-002	Barrel Clamp (Bushing)
326	321022	Barrel Clamp (Compression Spring)
158	TC10004194	Barrel Clamp (Park Position Collar) (Need to order Screw See Below)
N/A	10014496	Barrel Clamp (Park Position Collar NYLOK Set Screw)
335	147810-001	Barrel Clamp (Potentiometer Extension Guide)
321	302513	Barrel Clamp (Retaining Ring)
337	321155	Barrel Clamp (Shaft Seal)
340	147845-002	Barrel Clamp (Shaft Support)
365	TC10003525	Board, Display Assembly
360	TC10005079	Board, Flange Detector Sensor Assembly
395	147816-105	Board, IUI Interface Assembly
15	TC10008848	Board, Linear Sensor Assembly
390	TC10007649	Board, Logic Assembly (PCA Only) Use up TC10005342
8100	10010692	Calibration Kit, Syringe/PCA (Complete Cal Kit)
8100	10010691	Calibration Kit, Syringe/PCA (Force Sensor Spring Tool Only)
8100	148182-100	Calibration Kit, Syringe/PCA (Sizers and Height Gauge Tools)
9	49000214	Case, Front Assembly Kit (PCA Only)
406	147872-001	Case, Gasket Rear
8	10013672	Case, Rear Assembly Kit (PCA Only)
834	TC10006455	Door, Front (PCA Only)
832	148285-000	Door, Front Lockplate (PCA Only)
837	TC10003086	Door, Front Lockplate Strikeplate (PCA Only)
836	TC10003331	Door, Front Strikeplate (PCA Only)
	TC10003130	Door, Hinge (Bushings) (PCA Only)
835	TC10003128	Door, Hinge (PCA Only) (7 Screws #10012051 Sec Torx) (Lower Section of hinge requires 4 Hinge Bushings #TC10003130 & 4 Screws #10012051 Sec Torx for attaching to Front Case)
831	148287-000	Door, Rear (PCA Only)
838	148286-000	Door, Rear LockPlate (PCA Only)
350	TC10007202	Flag Barrel Position
355	147848-001	Flag Leaf Spring
345	TC10007599	Flag Seal
246	TC10003157	Flex Cable Retainer
250	10013882	Force Sensor Assembly Kit (includes screws NYLOCK #10013293)
203	147815-000	Force Sensor Seal

10	147539-101	Frame Drive Train Extrusion (3 Screws #321033 & 5 Screws #300923)
380	147739-100	Frame, Internal (2 Screws #301708 & 3 Screws #300923)
120	147537-000	Guide Rod
450	TC10003262	Handle, Left (PCA Only)
455	TC10003263	Handle, Right (PCA Only)
809	10013795	Handset, With Label (PCA Only) (Only included Alaris logo label 148349-000)
814	10210838	Handset, label "Patient Use Only" (Red label a like stop sign)
N/A	10942004	PCA Handset Cord Clip
804	148299-100	Harness, Handset Assembly Ground Strap
385	TC10008309	IUI Bracket
375	147792-000	IUI Bracket Strap (Ground Strap)
400	TC10008308	IUI Support Bracket
860	10010964	Key Kit (2 Keys per Kit) (PCA Only)
N/A	10014678	Kit, optional Locking Pole Clamp (with Locking knob Assy) Listed in DFU
373	147851-101	Keypad LED, Backlight Harness
2004	147078-100	Kit, IUI Connector Left
2003	147077-100	Kit, IUI Connector Right
2006	10013532	Kit, Latch (Consists of: Latch, Latch Support, Leaf Spring, Compression Spring, 2 Screws #10010865 Sec.Torx and Latch Security Bracket) (PCA Only)
5	49000225	Kit, Lead Screw Assembly
4	49000213	Kit, Lower Housing/Carriage assy. RoHS
645	148013-000	Label Patent
421	TC10003259	Latch, Bracket Module (PCA Only) (Security Bottom Latch Cover)
416	148470-000	Latch, Lever Module (PCA Only) (Security Door)
417	148494-000	Latch, Lever Module Cover (PCA Only) (Security Door)
800	148242-001	Latch, Module Bushing (PCA Only) (Latch Module Bottom Plate)
802	TC10003090	Latch, Module Bushing Support (PCA Only) (Latch Module Bottom Plate)
801	148245-001	Latch, Module Rod (PCA Only) (Latch Module Bottom Plate)
418	148495-000	Latch, Seal Module (PCA Only) (Security Door)
295	147741-000	Lens Status Indicator
205	147814-000	Lock Plate, Force Sensor
2001	10013524	Lock, Assembly Kit (PCA Only)
2	49000213	Housing Assy. Complete Assembly (Includes Force Sensor Board and Complete Actuator Knob Assembly) (Swap & Go Assembly)
6	49000248	Motor Pulley Assembly
85	TC10008310	Motor, (Bushing Isolation)
93	321050	Motor, (Rubber Grommet, Used with Motor Bushing)
461	321100	O-Ring, (Lower Housing to Upper Housing)
70	TC10008307	Plate Bottom Frame
200	TC10003887	Plate, Force Sensor
398	TC10003920	Retainer Clip for (Cable Assy, Dual End Connector)
410	142794-000	Rubber, Feet

480	303655	Screw, (Board, Flange Detector Sensor Assembly)
486	300923	Screw, (Board, Logic Assembly) & (Internal Frame at Bottom) & (Drive Train Assy. To Front Case) & (Latch Module Bracket)
498	10012051	Screw, (Door Hinge Security Torx) (PCA Only)
497	10012053	Screw, (Door, Front Lockplate Strikeplate) (PCA Only)
460	321033	Screw, (Drive Train Extrusion Frame Bottom Plate)
474	10015818	Screw, (Force Sensor NYLOCK)
488	301708	Screw, (IUI Bracket) & (Internal Frame at Top)
490	320855	Screw, (IUI Connector) & (Rear Case)
422	10010865	Screw, (Latch Module) Security Torx
478	148225-000	Screw, (Lower Housing to Upper Housing) & (Lock Plate, Force Sensor)
462	321060	Screw, (Motor Assembly)
492	320851	Screw, (Rear Case, Bottom Left) & (Latch Kit) & (Left Handle)
500	12276462	Screw, Dose request cord wrap
496	321035	Screw, (Right Handle)
327	12278652	Sensor Assembly, Syringe Size
396	806112	Snap Rivot
806	148373-100	Strap Wrap, Dose Request Assembly (Wrap & Brass Snap Will need to order screw #801131 and washer #320761) (PCA Only)
	12232252	Technical Service Manual (Syringe/PCA)
112	321031	Timing Belt
281	147093-100	Tubing MED Silicon, 10 feet
285	TC10007851	Upper Housing (4 Screws #321062 & 4 O-Rings #321100)
491	320729	Washer, (Rear Case)
494	320167	Washer, (Rear Case, Bottom Left)
808	305316	Washer, (#4 Internal Tooth Locking) Dose request cord