

525	142815-000	AIL (Gasket)
2001	49000295	AIL Sensor Kit, 8100
2007	49000204	Bezel Assembly (Mechanism Assembly) (Consists of: Bezel, O-Ring, Retaining Ring, Latch Roller, Compression Spring, Latch Yoke, Pin Yoke, Primary Finger, Secondary Finger, Occluder Fingers, Compression Springs, Camshaft Assembly, Mechanism Frame, Screws (7PL #320851), Ball Bearings, Cluster Gear, Torque Finger, Compression Springs, Torque Housing, and Gear Box) Refer to SB573 for old bezel kit 10964559 & door harness assy TC10006228 compatibility with the new bezel kit 49000204 & new door harness assy TC10007928
370	TC10006460	Bezel Assy Door Cover
220	TC10007928	Bezel Assy Door Harness (LVP Door Harness) Note: Old door harness TC10006228 is compatible with old bezel kit 10964559. New door harness TC10007928 is compatible with new bezel kit 49000204 (Refer to SB573)
442	TC10006587	Bezel Assy Membrane Frame
412	TC10004754	Board Display
470	12281195	Board Logic, Assy Kit (V9.1.17 LVP) Customer can use it with V8 LVP, but make sure they follow the insert card in the kit. Use up TC10006355
420	320129	Cable, Board Motor Controller & Logic Board Ribbon Cable
422	TC10008153	Board Motor Controller (New) Use up TC10005340 (CO 1077088)
423	148300-100	Board Motor Controller Kit (Motor & Motor Controller Board)
N/A	8100-PCS	Calibrated Set Pressure Control
N/A	8100-RCS	Calibrated Set Rate Control
2004	49000239	Case, Door Assy Kit (Front Door / Keypad Assy.)10360242 S/S by ROHS Kit per CO1098500
2009	10013325	Case, (Rear Assembly, Consists of: Nameplate Label, Patent Label, Rubber Foot/2, Kapton Tape, Silicon Tubing and Serial Number Replacement Label)
2003	49000222	Door Latch Assy. (Consists of: Door Latch, Pin, Sear, Setscrew, Latch Spring, and Sear Spring) (Now includes Spring & Latch Pivot Screw)
221	49000250	Door Screw (Latch Pivot)
451	320907	Door Spring
18	142794-000	Foot Rubber
2000	147082-100	General Hardware Kit, (Pump Module) (Consists of: Rubber Feet, Cluster Gear Ball Bearing, Spring Leaf Latch Module, Compression Spring, Snap Rivet, Washer #5 Flat, Washer Syringe Holder, Screw Latch, Screw Bezel Assy., Screw Rear Case, Screw IUI, Screw AIL Assv)
444	146332-100	Harness, Assembly, Pressure Transducer
220	TC10007928	Harness, Door Assembly. Use up TC10006228
410	TC10008302	IUI Bracket
2004	147078-100	IUI Connector Kit Left (Consists of: IUI Connector Left and Gasket)
2107	147077-100	IUI Connector Kit Right (Consists of: IUI Connector Right and Gasket)
643	10012610	Label, Patent
646	125569	Label, Serial Number Replacement
2006	147080-100	Latch Kit (Consists of: Latch, Latch Support, Leaf Spring and Compression Spring)
N/A	10848481	Manual Technical
444	146332-100	Mechanism Frame (Pressure Transducer Harness)
N/A	320079	Motor Controlled Board (Fuse F1) (.75a 125V)
2015	148438-100	Platen Assembly (Consists of: Hinge Pin, Hinge, and Torsion Spring)
750	143715-100	Pressure sensor Upper and Lower
482	146465-001	Pump Mech (Motor/Harness) (New)
482	146465-000	Pump Mech (Motor/Harness) (Old)

280	TC10006329	Pump Mech Assy Moutor Mount
261	12280874	Pump Mech Assy O-Ring
2005	147087-100	Pump Mech Assy Pinion Gear & Coupling Assembly (Gear & Coupling)
423	148035-001	Retainer Clip, (50 pin Dual Assembly Ribbon Cable) (Blue Clip)
443	321112	Retainer Clip, 5 Pin ALL Connector (Black Cap)
832	320848	Screw, (ALL)
815	320849	Screw, (Bezel Assy)
831	304319	Screw, (Board Motor Controller & Logic Board Screws)
830	320855	Screw, (IUI Connector) & (Rear Case)
800	320851	Screw, (Pump Mech Assy Gear Box Screw) & (Bottom of Latch)
825	320850	Screw, (Rear Case) & (Bezel Assy Door Cover)
N/A	10010865	Screw, (Security Torx used for secure locking of Module to PC Unit)
2009	147093-100	Silicone tubing (AKA: seal)
414	806112	Snap Rivot
	TC10005890	Status Indicator Lens
421	320789	Standoff, (ALL)
833	320908	Washer, Rear Case (Bottom)
445	148052-000	Washer, Rear Case (Top) (Washer Syringe Holder)
	8101A	Syringe Holder for LVP8100