

35	TC10005026	Board, (Display Masimo/Nellcor 8220/8210)
55	TC10007943	Board, (Logic Masimo/Nellcor 8220/8210) Use up 146637-106
37	10013006	Board, (Main Assembly Masimo 8220)
37	10014345	Board, (Main Assembly Nellcor 8210)
65	147119-100	Board, (Patient Connector Masimo 8220)
65	147767-100	Board, (Patient Connector Nellcor 8210)
39	146970-102	Board, (Power Supply Masimo 8220)
39	147747-100	Board, (Power Supply Nellcor 8210)
705	320930	Cable, (Patient, 3 Meter Nellcor)
705	320886	Cable, (Patient, 8 Ft Masimo)
20	TC10006216	Case, (Front Masimo 8220)
20	TC10006215	Case, (Front Nellcor 8210)
85	10014804	Case, Massimo Rear Case Assy Kit (SPO2 8220)
85	10014803	Case, Nellcor Rear Case Assy Kit (SPO2 8210)
75	142794-000	Foot Rubber
47	320906	Frame Assy Christmas Tree (Clips) 8220/8210
40	TC10008320	IUI Connector Bracket 8220/8210
60	146992-001	IUI Connector Bracket Strap 8220/8210
130	147078-100	IUI Connector Kit Left (IUI Connector Left and Gasket) 8220/8210
120	147077-100	IUI Connector Kit Right (IUI Connector Right and Gasket) 8220/8210
45	146977-000	IUI Frame 8220/8210
30	TC10005939	Keypad Assembly 8210/8220 (Carefusion logo, start using 9/1/2011)
32	144424-100	Keypad Backlight LED Assembly 8220/8210
631	147868-000	Label Defibrillation Rating 8210/8220
630	147381-000	Label Logo Masimo 8220
630	147755-000	Label Logo Nellcor 8210
620	10014353	Label Patent 8210
620	147352-000	Label Patent 8220
	147080-100	Latch Kit (Latch, Latch Support, Leaf Spring and Compression Spring)
	10015423	Manual Technical 8220/8210
36	320905	Nylon Standoff, (Board Power Supply/Masimo/Nellcor 8220/8210)
110	320871	Nylon Washer, (Board Power Supply/Masimo/Nellcor 8220/8210)
800	320851	Screw, (Board Display)
835	320855	Screw, (Board Power Supply/Masimo/Nellcor 8220/8210) & (IUI Connector)
825	300348	Screw, (Frame Assy Front Case) 8220/8210
830	806112	Snap Rivot
50	TC10006640	Speaker 8220/8210 (use 141117 until phased out)
15	146980-000	Status Indicator Gasket
25	146978-001	Status Indicator Lens