



The goal is simple:
the best possible outcome

BD™ Quincke Spinal Needles












BD™ Quincke spinal needles

BD has produced regional anesthesia products for most of our 118-year history. We have the experience, expertise and technology necessary to help you in your practice.



- Translucent hub provides contact clarity, which helps visualize cerebrospinal fluid (CSF) to ensure proper needle placement
- Key/slot arrangement of stylet and cannula hubs facilitates proper needle bevel orientation
- Sterile, single use

BD™ Quincke spinal needles			
Cat. no.	Description	Color	Packaging
405125	18G x 3.5" (90 mm)		25 / box
405123	20G x 3.5" (90 mm)		25 / box
405121	22G x 3.5" (90 mm)		25 / box
405120	23G x 3.5" (90 mm)		25 / box
405122	25G x 3.5" (90 mm)		25 / box
405124	26G x 3.5" (90 mm)		25 / box
405127	27G x 3.5" (90 mm)		25 / box
405126	27G x 3.5" (90 mm) - Set with Introducer		25 / box
405245	25G x 2" (51 mm) - Neonatal/Pediatric		25 / box



Disclaimer: The information provided herein is not meant to be used, nor should it be used, to diagnose or treat any medical condition. All content, including text, graphics, images and information etc., contained in or available through this literature is for general information purposes only. For diagnosis or treatment of any medical condition, please consult your physician/doctor. Becton Dickinson India Private Limited and or its affiliates, its employees are not liable for any damages/claims to any person in any manner whatsoever.

BD Medical

Becton Dickinson India Private Limited
5th & 6th Floor, Signature Tower – B
South City I, Gurugram, 122001
Haryana, India
Tel: 91-124-3949390
E-mail: bd_india@bd.com

bd.com

© 2018 BD. BD and the BD Logo are trademarks of Becton, Dickinson and Company.

BD-8458

