



# BD® Plastic NRFit™ Syringes

Trusted products, now available  
with ISO 80369-6 compliant connection

Misconnections can have devastating consequences – from permanent nerve injury<sup>1</sup> to paraplegia,<sup>2</sup> and even death.<sup>3,4</sup> Just one of these errors is simply one too many. That's why the International Standards Organization developed ISO 80369-6, a new design standard that requires neuraxial device connector tips to be incompatible with small-bore connectors used in other applications.

For over 100 years, BD has built its legacy on performance and trust through the delivery of high-quality products, service and support – resulting in a portfolio of hypodermic syringes that is more widely used by healthcare professionals globally than any other brand.

BD now offers a comprehensive range of NRFit™ ISO 80369-6 compliant syringes, available with slip and lock tips, for your regional anaesthesia and other neuraxial procedural needs.

BD® Plastic NRFit™ Syringes	
Reference	Description
400050	BD® 3 mL Syringe NRFit™ Lok
400172	BD® 3 mL Syringe NRFit™ Slip
400173	BD® 5 mL Syringe NRFit™ Lok
400051	BD® 5 mL Syringe NRFit™ Slip
400175	BD® 10 mL Syringe NRFit™ Slip (available mid-2021)
400174	BD® 10 mL Syringe NRFit™ Lok (available mid-2021)
400182	BD® 20 mL Syringe NRFit™ Lok
400183	BD® 50 mL Syringe NRFit™ Lok

For BD® Syringes NRFit™ Lok compatibility with infusion pumps, please refer to the pumps manufacturer's documentation.



To learn more, contact your BD representative

#### References

**1** Colak A, Ege T, Arar C, Yuksel V, Duran E. A severe complication of accidental epidural administration of glutaraldehyde. *Trakya Universitesi Tip Fakultesi Dergisi*. 2009;26(1):83–86. **2** Ajmal M. Accidental intrathecal injection of aminophylline in spinal anaesthesia. *Anesthesiology*. 2011;114(4):998–1000. **3** Koczmara C, Hyland S, Cheng R. Epidural medications given intravenously may result in death. *Dynamics*. 2007;18(3):34–36. **4** Lanigan, CJ. Safer epidural and spinal connectors. *Anaesthesia*. 2002;57(6):567–571



BD Switzerland Sàrl, Terre-Bonne Park - A4, route de Crassier 17, 1262 Eysins, Switzerland

[bd.com](http://bd.com)

