

# CME products accessories



BodyGuard  
**121**  
Twin™



BodyGuard  
**323**  
ColorVision™



BodyGuard  
**595**  
ColorVision™



**T34L**





## Accessories for




- |                        |      |
|------------------------|------|
| BodyGuard ColorVision™ | P. 2 |
| BodyGuard 121 Twins™   | P. 7 |
| T34, T34L, T34L-PCA    | P. 8 |



**CareFusion**

has joined BD





Item	Part Number	Description
<p>Rechargeable Li-Polymer Battery, 1800mAh</p> 	<p><b>150-050XV</b></p>	<ul style="list-style-type: none"> <li>• Standard battery supplied with BodyGuard ColorVision™.</li> <li>• Provides up to 17 hours at 125ml/h.</li> </ul>
<p>Rechargeable Li-Polymer Battery, 3600mAh</p> 	<p><b>130-051XV</b></p>	<ul style="list-style-type: none"> <li>• Optional "Extended life" Battery.</li> <li>• Provides up to 34 hours operating time at 125ml/h.</li> </ul>
<p>Power source 9 Volts x 2</p> 	<p><b>130-050X2V</b></p>	<ul style="list-style-type: none"> <li>• Holds (2) 9 Volt alkaline batteries.</li> <li>• Provides up to 12 hours pump operation at 125ml/h rate.</li> </ul>
<p>Charger BodyGuard ColorVision™</p> 	<p><b>150-320XV</b></p>	<ul style="list-style-type: none"> <li>• Docking station for BodyGuard ColorVision™ pumps.</li> <li>• Includes remote display, communication and drop sensor ports.</li> <li>• Works with a standard 230V power cord.</li> </ul>

Item	Part Number	Description
<p>Table top Charger</p> 	<p><b>150-331X</b></p>	<ul style="list-style-type: none"> <li>• Recharges standard and “extended life” Li-Polymer batteries.</li> <li>• Connects to wall (AC) outlet using standard power cord.</li> <li>• Battery recharge time approximately 6 to 8 hours for full capacity.</li> </ul>
<p>Wall charger European 8,2Volts</p> 	<p><b>151-142X</b></p>	<ul style="list-style-type: none"> <li>• Connects to wall (AC) power.</li> <li>• Connector plugs into DC Port on side of the BodyGuard ColorVision™ pump.</li> </ul>
<p>Patient White Handset with light</p> 	<p><b>140-400X</b></p>	<ul style="list-style-type: none"> <li>• Connects to Bolus port on the pump side.</li> <li>• Allow Patient bolus delivery.</li> </ul>

Item	Part Number	Description
<p>Lock Box 100ml</p> 	<p><b>170-100X</b></p>	<ul style="list-style-type: none"> <li>• Attaches to back of all Bodyguard single channel pumps.</li> <li>• Holds up to 100ml fluid bag.</li> <li>• Includes 2 Keys.</li> </ul>
<p>Lock Box 250ml</p> 	<p><b>170-250X</b></p>	<ul style="list-style-type: none"> <li>• Attaches to back of all Bodyguard single channel pumps.</li> <li>• Holds up to 250ml fluid bag.</li> <li>• Includes 2 Keys.</li> </ul>
<p>Lock Box ultimate 500ml</p> 	<p><b>190-900PXCVI</b></p>	<ul style="list-style-type: none"> <li>• Pole- mountable.</li> <li>• Holds up to 500ml fluid bag and a 50/60ml syringe.</li> <li>• Can fit with pump and charger and be placed on bed or pole.</li> <li>• Two locks: 1 for the pump/bag access and 1 to lock onto IV pole.</li> </ul>




Item	Part Number	Description
<p>Keys for 100ml and 250ml lock boxes</p> 	<p><b>170-010X</b></p>	<ul style="list-style-type: none"> <li>Compatible with lock box 100ml (170-100X) and lock box 250ml (170-250X).</li> </ul>
<p>Keys for 500ml ultimate lock box</p> 	<p><b>190-960XCV</b></p>	<ul style="list-style-type: none"> <li>Compatible with lock box ultimate 500ml (190-900PXCVI).</li> </ul>
<p>Pole Clamp adaptor</p> 	<p><b>CNC00002</b></p>	<ul style="list-style-type: none"> <li>Compatible with all BodyGuard, BodyGuard ColorVision™ pumps.</li> </ul>
<p>100ml Carrying bag blue</p> 	<p><b>150-000X</b></p>	<ul style="list-style-type: none"> <li>Holds BodyGuard ColorVision™ pump up to 100ml bag fluid volume.</li> </ul>




Item	Part Number	Description
<p>250ml Carrying bag</p> 	<p><b>Yellow: 150-250XE</b> <b>Blue: 100-083A</b></p>	<ul style="list-style-type: none"> <li>• Holds BodyGuard ColorVision™ pump up to 250ml bag fluid volume with or without lock box.</li> </ul>
<p>500ml Carrying bag</p> 	<p><b>150-001X</b></p>	<ul style="list-style-type: none"> <li>• Holds BodyGuard ColorVision™ up to 500ml bag (soft) and 250ml bag, semi-rigid.</li> </ul>
<p>2L TPN Carrying bag</p> 	<p><b>Grey: 160-002XG</b> <b>Black: 160-002XB</b></p>	<ul style="list-style-type: none"> <li>• Holds BodyGuard ColorVision™ pump up to 2l bag TPN solution.</li> <li>• Convenient for Paediatric patients.</li> </ul>
<p>3L TPN Carrying bag blue</p> 	<p><b>160-001X</b></p>	<ul style="list-style-type: none"> <li>• Holds BodyGuard ColorVision™ pump up to 3l bag TPN solution.</li> </ul>


Item	Part Number	Description
<p>Multi-dock</p> 	<p><b>100-100DSCV</b></p>	<ul style="list-style-type: none"> <li>• Holds 3 BodyGuard ColorVision™ pumps with their charger.</li> <li>• Fits on a pole.</li> </ul>
<p>Drop sensor</p> 	<p><b>121-100X</b></p>	<ul style="list-style-type: none"> <li>• Compatible with BodyGuard ColorVision and BG121™ Twins.</li> </ul>
<p>Cable Bodycomm USB-Serial</p> 	<p><b>197-100X</b></p>	<ul style="list-style-type: none"> <li>• Connects BodyGuard ColorVision™ charger with to a computer.</li> <li>• Includes CD for installation.</li> </ul>
<p>DC Power Cable BodyGuard</p> 	<p><b>196-000X</b></p>	<ul style="list-style-type: none"> <li>• Connects BodyGuard ColorVision™ to table top charger or pole mount charger while used in 100ml or 250ml lock box.</li> </ul>

Item	Part Number	Description
	<b>100-174S</b>	<ul style="list-style-type: none"> <li>• Fits with T34 syringe pump.</li> <li>• Delivered with 2 keys.</li> </ul>
	<b>100-174SL</b>	<ul style="list-style-type: none"> <li>• Fits with T34L syringe pump.</li> <li>• Delivered with 2 keys.</li> </ul>
	<b>100-176S</b>	<ul style="list-style-type: none"> <li>• Compatible with T34 syringe pump fitted with lock box (100-174S).</li> </ul>
	<b>100-178S</b>	<ul style="list-style-type: none"> <li>• Compatible with T34L syringe pump stand alone.</li> </ul>



Item	Part Number	Description
<p>Carrying bag T34 no Lock Box</p> 	<p><b>100-177S</b></p>	<ul style="list-style-type: none"> <li>• Compatible with T34 syringe pump stand alone.</li> </ul>
<p>Patient blue Handset with light</p> 	<p><b>140-100X</b></p>	<ul style="list-style-type: none"> <li>• Connects to Bolus port on side of the pump (T34L-PCA).</li> <li>• Allow Patient bolus delivery.</li> <li>• Compatible with BodyGuard ColorVision™ pumps.</li> </ul>
<p>Key for T34 lock box USB-Serial</p> 	<p><b>100-206SL</b></p>	<ul style="list-style-type: none"> <li>• Compatible with T34 lock box (100-174S).</li> </ul>

Item	Part Number	Description
<p>Key for lock box</p> 	<p><b>190-115X</b></p>	<ul style="list-style-type: none"> <li>• Compatible with T34L &amp; T34L-PCA lock box (100-1745L).</li> </ul>
<p>Docking station for T34L and T34L-PCA</p> 	<p><b>130-000DS</b></p>	<ul style="list-style-type: none"> <li>• Fits to a pole.</li> <li>• Connects to wall (AC) power plug.</li> <li>• Hold T34L and T34L-PCA pump.</li> <li>• Can be stacked vertically up to 7 pumps.</li> </ul>
<p>Wall charger European 8,2Volts</p> 	<p><b>151-142X &amp; 151-141X</b></p>	<ul style="list-style-type: none"> <li>• Connects to wall (AC) power</li> <li>• Connector plugs into DC Port on side of T34L and T34L-PCA</li> <li>• 151-141X (8.4V) Recommended for ambulatory use.</li> <li>• 151-142X (8,2V) Recommended for stationary use.</li> </ul>

Item	Part Number	Description
<p>Communication station T34</p> 	<p><b>150-312S</b></p>	<ul style="list-style-type: none"> <li>• Enables T34 syringe pump connection to a computer.</li> <li>• Battery operated.</li> </ul>

Manufactured by CME

**CE**  
0344

EC REP

**MedNet GmbH**  
Brokstrasse 10 48163  
Muenster, Germany  
+49 251 32266 0

**CME**

Caesarea Medical Electronics



**Caesarea Medical Electronics Ltd.**  
16 Shacham Street, Industrial Park  
Caesarea North, P.O.BOX 4294  
Caesarea 38900, Israel  
[www.cme-infusion.com](http://www.cme-infusion.com)

**CareFusion Switzerland 317 Sarl**  
A-One Business Centre, Z.A.  
Vers-La-Pièce n°10  
CH-1180  
Rolle

© 2015 CareFusion Corporation or one of its subsidiaries  
Distributed by CareFusion, Please contact your local representative. 0000CF02552 Issue1

[carefusion.com](http://carefusion.com)

Distributed by



**CareFusion**  
*has joined BD*