



Pyxis SupplyCenter server

Pyxis® Perioperative Solutions

The Pyxis® SupplyCenter server collects data from Pyxis supply management technologies located throughout your facility or IDN, allowing you to track and analyze patient, item, user and procedure-specific supply information per facility or across an IDN. Acting as the central data repository for all Pyxis supply automation activity, the Pyxis SupplyCenter server lets you keep an eye on your bottom line while helping ensure clinicians have what they need at the point of care.

Pyxis SupplyCenter optimizes supply management by:

- Consolidating all point-of-use transactions to one location
- Accurately tracking reorder and inventory needs
- Fulfilling replenishment orders automatically
- Providing immediate notification of stock-outs
- Sending electronic billing updates to the hospital's billing system
- Sending usage information to the hospital's MMIS and other integral systems
- Providing integration with ORIS and other clinical systems

Pyxis®



Enterprise inventory and data management

Pyxis SupplyCenter is the central hub for all Pyxis supply automation activity—providing you quick access to information necessary for efficient supply management and optimal decision making.

Key advantages

- Generates more than 100 targeted reports and analyzes data per facility, or throughout an entire healthcare network
- Provides quick and easy management of inventory stock-outs, and reorders items directly from the system
- Dynamic Order Point feature recommends par level adjustments for tighter inventory control and ongoing cost savings
- Retains 45 days of transaction history and 13 months of summary data for each station, including supply utilization, refill activity, inventory verification, patient charges and credits
- Request Item feature from the station serves as a communication tool between the nursing unit and distribution, saving valuable clinical time and minimizing urgent phone calls
- Non-stock Item feature can be used to manage special order items, eliminating the need for clinical management
- Inactive Patient feature allows users to search for patients by medical record number or based on last activity time and discharge, as appropriate



Pyxis SupplyCenter configuration: LCD panel, CPU, keyboard, laser printer and rack.

Dimensions: 49" H x 30" W x 24" D; 128 lbs.

Available with or without desk.

A rack mountable blade server is an available option.

Pyxis SupplyCenter can support more than 100 stations and can accommodate up to 10 client workstation PCs.

CareFusion
San Diego, CA

carefusion.com

© 2011 CareFusion Corporation or one of its subsidiaries. All rights reserved. Pyxis is a trademark or registered trademark of CareFusion Corporation or one of its subsidiaries. PS2214 (0411/2000)

