

V. Mueller® instruments Cosgrove flex clamps

Cosgrove flex clamps from V. Mueller instruments have flexible shafts that bend and contour, yet they perform with the reliability of traditional rigid clamps.

V. Mueller®



Two flex clamp choices

We offer two choices in Cosgrove flex clamps. When a tighter bending radius and maximum flexibility are important, choose the Quick Bend™ model. The Quick Bend flex clamp has a shortened neck that bends closer to the jaw for a tighter bending radius to give you even greater flexibility.

And, when the width of the clamp head is key for working in small incisions, the Slim Fit™ model is the right choice for your needs. The Slim Fit flex clamp has a slender clamp head that permits use in narrow incisions.

To activate, just position and ratchet. The flex clamp becomes rigid and locks in place.

Cosgrove flex clamp features

- Strong, flexible, reliable shaft constructed of beaded stainless steel
- Enough clamping force for total aorta occlusion without excessive pressure
- Wide jaw opening for a secure grasp of all vessels
- Gap between shaft beads to allow for thorough cleaning
- Familiar ratchet design handle
- Slim design excellent for use in minimally invasive procedures





The power of partnership

The Cosgrove flex clamp was developed in conjunction with Dr. Delos Cosgrove and the Cleveland Clinic through the V. Mueller Custom Signature Service® program. The Cosgrove flex clamp is just another example of our commitment to innovation and to working closely with surgeons.

V. Mueller instrumentation—the name you trust

Cosgrove flex clamps are instruments you can use with confidence because they come from the V. Mueller division of CareFusion, one of the world's leading providers of surgical instruments for more than a century. The full line of Cosgrove flex clamps is beautifully designed and precision engineered for quality, performance and reliability.

V. Mueller surgical instrument guarantee

We warrant every surgical instrument bearing the V. Mueller brand name to be free of functional defects in workmanship and materials when used normally for its intended surgical purposes. Any V. Mueller instrument proving to be so defective will be replaced or repaired at our option, at no charge.






To place an order or for more information, please contact your V. Mueller Sales Representative or call 800.323.9088.

Quick Bend flex clamps	
Part no.	Description
CV1133	Cosgrove flex clamp 33 mm, Fogarty-type jaw
CV1161	Cosgrove flex clamp 61 mm, Fogarty-type jaw
CV1186	Cosgrove flex clamp 86 mm, Fogarty-type jaw

Clamp inserts not included.

Slim Fit flex clamps	
Part no.	Description
CV1033	Cosgrove flex clamp 33 mm, Fogarty-type jaw
CV1061	Cosgrove flex clamp 61 mm, Fogarty-type jaw

Clamp inserts not included.



To place an order or for more information, please contact your V. Mueller Sales Representative or call 800.323.9088.

CareFusion Waukegan, IL

carefusion.com

