



## BD SafetyGlide™

Easy activation keeps healthcare workers protected and efficient.



# Committed to creating a safer environment for both you and your patients.


In today's world of health, increasing workload and productivity means ensuring safety is more challenging than ever.

BD SafetyGlide™ needle, extends the available safety technologies to cover your needs.



In one hospital study, the use of BD SafetyGlide™ syringes and needles – together with blunt fill needles, training and education – resulted in a significant decrease in reported needlestick injuries from 16.9 to 6 / 100,000 devices<sup>1</sup>.

# Designed for your safety and to fit your natural workflow.



**Fluid single-finger activation** requiring minimal training and virtually no change in clinician technique

**Stainless steel latch** securely shields needle after activation

**Colour coded hub** for easy needle size identification

**Unique Activation-Assist feature** meaning a fast and smooth activation while keeping your finger well behind the needle

**Same needle sharpness** as conventional

**Bevel up needle orientation** to support recommended injection technique

# Ordering information

## BD SafetyGlide™ Safety Needles

Gauge size	Outer Diameter (mm)	Length (mm)	Length (inch)	Color Code	Shelf box / Shipping case	Reference number
27 G	0.4	16	5/8"	Grey	50 / 500	305921
25 G	0.5	16	5/8"	Orange	50 / 500	305901
25 G	0.5	25	1"	Orange	50 / 500	305916
23 G	0.6	25	1"	Blue	50 / 500	305902
22 G	0.7	40	1 1/2"	Black	50 / 500	305900
21 G	0.8	25	1"	Green	50 / 500	305915
21 G	0.8	40	1 1/2"	Green	50 / 500	305917
18 G	1.2	40	1 1/2"	Pink	50 / 500	305918

#### References:

1. Adams D, Elliott TS. Impact of safety needles devices on occupationally acquired needlestick injuries: A four-year prospective study. *J Hosp Infect.* 2006. Sept;64(1):50-5.

BD, Terre-Bonne Park – A4, Route de Crassier 17, 1262 Eysins, Switzerland  
Tel: +41 21 556 30 00

[bd.com](http://bd.com)

