


EUCAST



Win with the First Complete Offer

New BD Phoenix™ EUCAST compliant panels



Helping all people
live healthy lives

BD Phoenix™ EUCAST compliant panels



BD is committed to continuously providing its customers with up to date solutions. EUCAST is already a reality in the Nordic countries, UK and France, and we believe that EUCAST will soon be adopted as THE European Gold Standard.



We are proud to announce the availability of the first BD Phoenix™ EUCAST compliant panel set.

Phoenix Gram negative

Cat. No: 448103

Format: NMIC/ID-76

Phoenix Gram Negative		
code	Antimicrobial	(µg/mL)
AN	Amikacin	4 - 16
AXC	Amoxicillin/Clavulanate	2/2 - 8/2
AM	Ampicillin	2 - 8
ATM	Aztreonam	1 - 16
CZ	Cefazolin	1 - 4
CTX	Cefotaxime	0.5 - 4
CAZ	Ceftazidime	0.5 - 8
CXM	Cefuroxime	1 - 8
CN	Cephalexin	2 - 16
CIP	Ciprofloxacin	0.125 - 1
CL	Colistin	1 - 4
ETP	Ertapenem	0.25 - 1
ESR	ESBL	yes
FF	Fosfomycin	16 - 64
GM	Gentamicin	1 - 4
IPM	Imipenem	1 - 8
MEM	Meropenem	1 - 8
FM	Nitrofurantoin	16 - 64
TZP	Piperacillin/Tazobactam	4/4 - 16/4
NN	Tobramycin	1 - 4
TMP	Trimethoprim	1 - 4
SXT	Trimethoprim/ Sulfamethoxazole	1/19 - 4/76

Phoenix Gram positive

Cat. No: 448089

Format: PMIC/ID-67

Phoenix Gram Positive		
code	Antimicrobial	(µg/mL)
AXC	Amoxicillin/Clav.	2/2 - 8/2
AM	Ampicillin	2 - 8
NCF	Beta-lactamase	yes
FOX	Cefoxitin	2 - 8
CIP	Ciprofloxacin	0.5 - 2
CC	Clindamycin	0.25 - 1
DAP	Daptomycin	1 - 4
E	Erythromycin	0.25 - 2
FA	Fusidic Acid	0.5 - 8
GM	Gentamicin	1 - 4
GMS	Gentamicin-Synergy	500
LZD	Linezolid	0.5 - 4
MXF	Moxifloxacin	0.125 - 1
MUP	Mupirocin	1 - 4
MUH	Mupirocin - High Level	256
FM	Nitrofurantoin	16 - 64
OX	Oxacillin	0.25 - 2
P	Penicillin	0.0625 - 0.25
RA	Rifampin	0.25 - 1
TEC	Teicoplanin	1 - 8
TE	Tetracycline	0.5 - 2
NN	Tobramycin	1 - 4
TMP	Trimethoprim	1 - 4
SXT	Trimethoprim/ Sulfamethoxazole	1 / 19 - 4/76
VA	Vancomycin	0.5 - 8

Please contact your local BD representative for more details.



BD Diagnostics

11, rue Aristide Bergès
38800 Le Pont de Claix,
France

Tel: (33) 476 68 36 36

www.bd.com/diagnostics