


IMAGINE the open possibilities



BD MAX[™]

The Power of Choice





imagine the possibilities

BD MAX™ automates and standardizes your work flow

- Prepare samples
- Create onscreen work list
- Load rack(s)
 - Samples
 - Unitized reagent strips
 - User-defined primers/probes
- Place rack(s) in BD MAX[™]
- Load cartridge(s) and close door



Add user-defined primers/probes to unitized reagent strips — load rack(s)



Place rack(s) in BD MAX™



Load cartridge(s) and close door

BD MAX™

OPEN SYSTEM REAGENTS AND CONSUMABLES



BD MAX™ DNA Extraction Kits

437523 BD MAX DNA Extraction Kit DNA-1 (Urine/Plasma)

437524 BD MAX DNA Extraction Kit DNA-2 (Whole Blood)

437525 BD MAX DNA Extraction Kit DNA-3 (Swabs in Transport Medium/UTM)



BD MAX™ RNA Extraction Kits

437522 BD MAX RNA Extraction Kit RNA-3 (Swabs in Transport Medium/UTM)



BD MAX™ Generic Reagents

435014 BD MAX Master Mix MM (With Primers/Probe for Process Control)

435022 BD MAX Master Mix MM1

BD MAX™ Consumables

437519 BD MAX Microfluidic Cartridges, 24 Lane

For more information visit our website: www.bd.com/europe/ds/



BD Diagnostics Tullastrasse 8-12 69126 Heidelberg Germany Tel. (49) 6221 305 0 Fax. (49) 6221 305 604 www.bd.com/europe/ds/