

BD Beaver™ and BD Visitec™

Specialty Surgical Products



Helping all people
live healthy lives

- Ear, Nose, and Throat (ENT)
- Emergency Medical
- Neuro
- Orthopedic
- Other Specialty Surgical

Introduction

BD is committed to improving the quality of healthcare for all people worldwide. With a long-standing tradition of innovation in healthcare, BD is devoted to developing high quality products that meet a broad spectrum of surgical and clinical needs.

Our products are manufactured in the United States and United Kingdom, using the highest levels of care and quality. BD Beaver™ specialty knives, blades, and accessories are designed for Ear, Nose and Throat (ENT), Emergency Medical (EMT), Gynecology, Neuro Surgery, Orthopedic, and Reconstructive procedures.

Recognized throughout the industry for maintaining a commitment to quality, and safety, BD embraces the commitment by delivering products that have been designed to consistently meet or exceed industry standards. BD provides products, services, and educational programs to professionals, and their patients, in an effort to raise the standard of care for all the world.

At BD, our employees have the knowledge and technology to build partnerships with surgical professionals; creating solutions for today's surgical challenges. BD is making the necessary investments to ensure that we can continue meeting future challenges.

We appreciate your support and comments to help us all achieve better healthcare.

Visit our website at www.bd.com to see what's new and upcoming at BD.

www.bd.com

BD Beaver™ and BD Visitec™ Products

Table of Contents

EAR, NOSE, AND THROAT (ENT)

Mini-Blades and Knife	2
Myringotomy Blades and Knives	3
Needle Blades.	2
Tonsillectomy Blade	2
Tympanoplasty and Stapedectomy Blades	3

EMERGENCY MEDICAL (EMT)

Micro-Blades	5
Micro-Sharp™ Blades and Micro-Unitome™ Knives	5
Mini-Blades and Knife	4
Ring Cutter and Blades	4

NEURO SURGERY

Neuro-Sharp™ Blades	6
Needle and Sickle Edge Blades	6

ORTHOPEDIC

Open Procedures

Debakey Blades	8
Orthopedic Blades	7, 8

Arthroscopic Procedures

Arthro-Lok™ Blades	9
------------------------------	---

OTHER SPECIALTY SURGICAL PRODUCTS

Biopsy Blade	12
Debakey and Mammatome Blades	12
Disposable Instruments	12
Micro-Blades	11
Micro-Sharp™ Blades and Micro-Unitome™ Knives	11
Mini-Blades and Knife	10
Needle Blades	10

BD BEAVER™ HANDLES

Handles	13
-------------------	----

HANDLE INDEX	14, 15
------------------------	--------

PRODUCT INDEX	16, 17
-------------------------	--------

Ear, Nose, Throat (ENT) Products

BD BEAVER™ MINI-BLADES AND KNIFE

376300 Mini-Blade

sharp on two sides, pointed tip
Quantity - 6 per box

**376400 Mini-Blade**

rounded tip, sharp on one side
Quantity - 24 per box

**376500 Mini-Blade**

smaller alternative to #11 general blade
Quantity - 24 per box

**376700 Mini-Blade**

smaller alternative to #15 general blade
Quantity - 24 per box

**376900 Mini-Blade**

straight, sharp all around, 180° cutting surface
double bevel blade
Quantity - 6 per box

**376910 Mini-Blade**

angled 60°, sharp all around, 180° cutting surface
double bevel blade
Quantity - 6 per box

**374769 Mini-Blade Knife**

straight, sharp all around, 180° cutting surface
double bevel blade
Quantity - 3 per box



BD BEAVER™ NEEDLE BLADES

375910 Needle Blade

3.0 mm cutting edge
Quantity - 3 per box

**375920 Needle Blade**

5.0 mm cutting edge
Quantity - 3 per box

**375930 Needle Blade**

7.0 mm cutting edge
Quantity - 3 per box

**375940 Needle Blade**

1.5 mm cutting edge
Quantity - 3 per box



BD BEAVER™ OROPHARYNGEAL BLADES

376320 Interdental Blade

60° angled, sharp all around, pointed tip
Quantity - 6 per box

**376800 Tonsillectomy Blade**

sickle sharp
Quantity - 24 per box



Ear, Nose, Throat (ENT) Products

BD BEAVER™ MYRINGOTOMY BLADES AND KNIVES

377100 Myringotomy Blade

lance tip, downward cut, offset 45°
Quantity - 6 per box

**377110 Myringotomy Blade**

lance tip, upward cut, offset 45°
Quantity - 6 per box

**377120 Myringotomy Blade**

spear tip, narrow shaft, offset 45°
Quantity - 6 per box

**377121 Myringotomy Blade**

lance tip, downward cut, narrow shaft, offset 45°
Quantity - 6 per box

**377122 Myringotomy Knife**

spear tip, narrow shaft, offset 45°
Quantity - 3 per box

**377124 Myringotomy Knife**

lance tip, downward cut, narrow shaft, offset 45°
Quantity - 3 per box



BD BEAVER™ TYMPANOPLASTY AND STAPEDECTOMY BLADES

374451 Tympanoplasty Blade

50° angled, sharp all around
2.0 mm W, bevel down
Quantity - 3 per box

**377200 Tympanoplasty Blade**

60° angled, sharp all around
2.5 mm W, bevel down
Quantity - 3 per box

**377210 Tympanoplasty Blade**

straight, sharp all around
2.5 mm W
Quantity - 3 per box

**377300 Sickle Edge Blade**

straight 84 mm L
5.0 mm cutting surface
Quantity - 3 per box



Emergency Medical Products

BD BEAVER™ RING CUTTER AND BLADES

374100 Ring cutter

Quantity - 1 per box



374150 Ring cutter saw blade

Quantity - 1 per box



374151 Ring cutter saw blade

high speed steel

Quantity - 1 per box



BD BEAVER™ MINI-BLADES AND KNIFE

374769 Mini-Blade Knife

straight, sharp all around, 180°
cutting surface, double bevel blade

Quantity - 3 per box



376300 Mini-Blade

sharp on two sides, pointed tip

Quantity - 6 per box



376400 Mini-Blade

sharp on one side, rounded tip

Quantity - 24 per box



376500 Mini-Blade

smaller alternative to #11 general blade

Quantity - 24 per box



376600 Mini-Blade Lamellar Blade

angled 60°, smooth undersurface, sharp all around, bevel up

Quantity - 3 per box



376610 Mini-Blade Lamellar Blade

straight, smooth undersurface, sharp all around

Quantity - 3 per box



376700 Mini-Blade

smaller alternative to #15 general blade

Quantity - 24 per box



376900 Mini-Blade

straight, sharp all around, 180° cutting surface, double bevel blade

Quantity - 6 per box



Emergency Medical Products

BD BEAVER™ MICRO-SHARP™ BLADES AND MICRO-UNITOME™ KNIVES

377511 Micro-Sharp 15°
1.5 mm depth, orange
Quantity - 6 per box



377512 Micro-Unitome 15°
1.5 mm depth, orange
Quantity - 10 per box



377513 Micro-Sharp 15°
3.0 mm depth, blue
Quantity - 6 per box



377514 Micro-Unitome 15°
3.0 mm depth, blue
Quantity - 10 per box



377515 Micro-Sharp 15°
5.0 mm depth, green
Quantity - 6 per box



377516 Micro-Unitome 15°
5.0 mm depth, green
Quantity - 10 per box



377522 Micro-Sharp 22.5°
4.0 mm depth, light blue
Quantity - 6 per box



377523 Micro-Unitome 22.5°
4.0 mm depth, light blue
Quantity - 10 per box



377530 Micro-Sharp 30°
3.5 mm depth, black
Quantity - 6 per box



377531 Micro-Unitome 30°
3.5 mm depth, black
Quantity - 10 per box



377545 Micro-Sharp 45°
3.0 mm depth, tan
Quantity - 6 per box



377546 BD Micro-Unitome 45°
3.0 mm depth, tan
Quantity - 10 per box



BD BEAVER™ MICRO-BLADES

377615 Micro-Blade
Angled 15°
Quantity - 25 per box



377630 Micro-Blade
Angled 30°
Quantity - 25 per box



Neuro Surgery Products

BD BEAVER™ NEURO-SHARP™ BLADES

370110 Neuro-Sharp Blade

15° offset sickle

Quantity - 3 per box



370120 Neuro-Sharp Blade

15° offset arachnoid

Quantity - 3 per box



370410 Neuro-Sharp Blade

bayonet sickle

Quantity - 3 per box



370420 Neuro-Sharp Blade

bayonet arachnoid

Quantity - 3 per box



BD BEAVER™ NEEDLE AND SICKLE EDGE BLADES

375910 Needle Blade

3.0 mm cutting edge

Quantity - 3 per box



375920 Needle Blade

5.0 mm cutting edge

Quantity - 3 per box



375930 Needle Blade

7.0 mm cutting edge

Quantity - 3 per box



375940 Needle Blade

1.5 mm cutting edge

Quantity - 3 per box



377300 Sickle Edge

straight 84 mm L

5.0 mm cutting surface

Quantity - 3 per box



BD BEAVER™ ORTHOPEDIC BLADES

371400 Bent Blade

sharp all around, bent 45°
Quantity - 4 per box

**371600 Bent Blade**

sharp one edge, bent 45°
Quantity - 4 per box

**371800 Bent Blade**

sharp two edges, blunt end, bent 45°
Quantity - 4 per box

**374562 Mini-Blade**

Quantity - 24 per box

**376100 Mini-Blade**

Quantity - 24 per box

**378100 Chisel Blade**

non-sterile
Quantity - 5 per box

**378300 Chisel Blade**

non-sterile
Quantity - 5 per box

**378400 Chisel Blade**

non-sterile
Quantity - 5 per box

**378600 Chisel Blade**

non-sterile
Quantity - 5 per box

**378800 Chisel Blade**

non-sterile
Quantity - 5 per box



Orthopedic Open Procedure Products

BD BEAVER™ ORTHOPEDIC BLADES

379400 Chondroplastic Blade
sharp on two sides, curved
Quantity - 4 per box



379500 Meniscus Blade
curved, 7 mm
Quantity - 4 per box



379700 Meniscus Blade
curved, 9 mm
Quantity - 4 per box



379710 Meniscus Blade
flat, 9 mm
Quantity - 4 per box



BD BEAVER™ DEBAKEY BLADES

374410 Blade
5 mm cutting width
Quantity - 24 per box



374710 Blade
8 mm cutting width
Quantity - 24 per box



374810 Blade
10 mm cutting width
Quantity - 24 per box



Orthopedic Arthroscopic Procedure Products

BD BEAVER™ ARTHRO-LOK™ BLADES

374580 Arthro-Lok Blade

4 mm hoe, sharp all around, angled 60°
Quantity - 6 per box



376083 Arthro-Lok Blade

4 mm retrograde
Quantity - 6 per box



376086 Arthro-Lok Blade

4 mm pointed retrograde
Quantity - 6 per box



376480 Arthro-Lok Blade

4 mm rounded tip, curved left, sharp on one side
Quantity - 6 per box



376481 Arthro-Lok Blade

4 mm rounded tip, flat, sharp on one side
Quantity - 6 per box



376482 Arthro-Lok Blade

4 mm rounded tip, curved right, sharp on one side
Quantity - 6 per box



376780 Arthro-Lok Blade

4 mm pointed tip, curved left, sharp on one side
Quantity - 6 per box



376782 Arthro-Lok Blade

4 mm pointed tip, curved right, sharp on one side
Quantity - 6 per box



376980 Arthro-Lok Blade

4 mm mini-banana, sharp all around, curved
Quantity - 6 per box



376981 Arthro-Lok Blade

4 mm mini-banana, sharp all around, flat
Quantity - 6 per box



376984 Arthro-Lok Blade

4 mm banana, sharp all around
Quantity - 6 per box



377081 Arthro-Lok Blade

4 mm rosette
Quantity - 6 per box



379071 Arthro-Lok Blade

4 mm crescent, flat
Quantity - 6 per box



379081 Arthro-Lok Blade

4 mm mini-meniscus, flat
Quantity - 6 per box



379082 Arthro-Lok Blade

4 mm mini-meniscus, curved
Quantity - 6 per box



Other Specialty Surgical Products

BD BEAVER™ MINI-BLADES AND KNIFE

374769 Mini-Blade Knife

straight, sharp all around, 180°
cutting surface, double bevel blade
Quantity - 3 per box



376300 Mini-Blade

sharp on two sides, pointed tip
Quantity - 6 per box



376400 Mini-Blade

rounded tip, sharp on one side
Quantity - 24 per box



376500 Mini-Blade

smaller alternative to #11 general blade
Quantity - 24 per box



376600 Mini-Blade Lamellar Blade

angled 60°, smooth undersurface, sharp all around, bevel up
Quantity - 3 per box



376610 Mini-Blade Lamellar Blade

straight, smooth undersurface, sharp all around
Quantity - 3 per box



376700 Mini-Blade

smaller alternative to #15 general blade
Quantity - 24 per box



376900 Mini-Blade

straight, sharp all around, 180° cutting surface, double bevel blade
Quantity - 6 per box



376910 Mini-Blade

angled 60°, sharp all around, 180° cutting surface
double bevel blade
Quantity - 6 per box



BD BEAVER™ NEEDLE BLADES

375910 Needle Blade

3.0 mm cutting edge
Quantity - 3 per box



375920 Needle Blade

5.0 mm cutting edge
Quantity - 3 per box



375930 Needle Blade

7.0 mm cutting edge
Quantity - 3 per box



375940 Needle Blade

1.5 mm cutting edge
Quantity - 3 per box



Other Specialty Surgical Products

BD BEAVER™ MICRO-SHARP™ BLADES AND MICRO-UNITOME™ KNIVES

377511 Micro-Sharp 15°
1.5 mm depth, orange
Quantity - 6 per box



377512 Micro-Unitome 15°
1.5 mm depth, orange
Quantity - 10 per box



377513 Micro-Sharp 15°
3.0 mm depth, blue
Quantity - 6 per box



377514 Micro-Unitome 15°
3.0 mm depth, blue
Quantity - 10 per box



377515 Micro-Sharp 15°
5.0 mm depth, green
Quantity - 6 per box



377516 Micro-Unitome 15°
5.0 mm depth, green
Quantity - 10 per box



377522 Micro-Sharp 22.5°
4.0 mm depth, light blue
Quantity - 6 per box



377523 Micro-Unitome 22.5°
4.0 mm depth, light blue
Quantity - 10 per box



377530 Micro-Sharp 30°
3.5 mm depth, black
Quantity - 6 per box



377531 Micro-Unitome 30°
3.5 mm depth, black
Quantity - 10 per box



377545 Micro-Sharp 45°
3.0 mm depth, tan
Quantity - 6 per box



377546 Micro-Unitome 45°
3.0 mm depth, tan
Quantity - 10 per box



BD BEAVER™ MICRO-BLADES

377615 Micro-Blade
Angled 15°
Quantity - 25 per box



377630 Micro-Blade
Angled 30°
Quantity - 25 per box



Other Specialty Surgical Products

BD BEAVER™ DEBAKEY AND MAMMATOME BLADES

374410 Blade

5 mm cutting width
Quantity - 24 per box

**374710 Blade**

8 mm cutting width
Quantity - 24 per box

**374810 Blade**

10 mm cutting width
Quantity - 24 per box

**374900 Blade**

270° cutting edge
Quantity - 6 per box



BD BEAVER™ GYNECOLOGY BLADE

379100 Cervical Biopsy Blade

Quantity - 4 per box



BD VISITEC™ DISPOSABLE INSTRUMENTS

581409 Bonn Forceps (0.12 tip)

Disposable instrument used to grasp tissues.
10 cm long - 12 mm long tip. 10 per box.

**581411 Tying Forceps**

Disposable instrument used in tying sutures. 10 cm
bent tip with micro-trays designed for 10/0 sutures.
10 per box.

**581417 Needle Holder**

Disposable instrument used to hold the needle
during suturing. 120 mm. 10 per box.

**581422 Vannas Scissors**

Disposable instrument. Small, sharp scissors used to
perform delicate cutting. 120 mm. 10 per box.



BD Beaver™ Handles

BD BEAVER™ HANDLES

Stainless steel handles available in round and hex configurations in a variety of lengths. All handles will accommodate all BD blades and collets. All handles and collets are non-sterile.

374311 Handle

Round, Knurled, Stainless, 5.0 cm L x .50 cm D

374310 Handle

Round, Knurled, Stainless, 10.0 cm L x .50 cm D

371350 Handle

Round, Knurled, Stainless, 13.0 cm L x .50 cm D

371360 Handle

Round, Knurled, Stainless, 13.0 cm L x .66 cm D

371394 Handle

Round, 18 cm L x .4 cm

371395 Handle

Round, 21 cm L x .4 cm

374320 Handle

Round, Knurled, Stainless, 15.5 cm L x .50 cm D

371410 Handle

Round, Knurled, Stainless, 10.2 cm L x .64 cm D

374330 Handle

Hex, Stainless, 10.0 cm L x .50 cm D

374340 Handle

Hex, Stainless, 15.5 cm L x .50 cm D

372000 BeaverTail™ Handle

Flat, 17 cm L

375000 BeaverTail Handle


Flat, 13 cm L

371305 Replacement Chuck Collet

Stainless, fits .50 cm handles

371415 Replacement Chuck Collet

Stainless, fits .60 cm handles



All handles - 1 per box

www.bd.com

Handles and Accessories

REPLACEMENT COLLETS

HANDLES

371350 371394 371395 372000 374310 374311 374320 374330 374340 375000 371305

BLADES

	371350	371394	371395	372000	374310	374311	374320	374330	374340	375000	371305
370110	x				x	x	x	x	x		x
370120	x				x	x	x	x	x		x
370410	x				x	x	x	x	x		x
370420	x				x	x	x	x	x		x
371400				x							x
371600				x							x
371800				x							x
374410					x					x	x
374451	x				x	x	x	x	x	x	x
374562	x				x	x	x	x	x		x
374580		x	x								
374710				x						x	
374810				x						x	
374900				x						x	
375910	x				x	x	x	x	x		x
375920	x				x	x	x	x	x		x
375930	x				x	x	x	x	x		x
375940	x				x	x	x	x	x		x
376083		x	x								
376086		x	x								
376100	x				x	x	x	x	x		x
376300	x				x	x	x	x	x		x
376400	x				x	x	x	x	x		x
376480		x	x								
376481		x	x								
376482		x	x								
376500	x				x	x	x	x	x		x
376600	x				x	x	x	x	x		x
376610	x				x	x	x	x	x		x
376700	x				x	x	x	x	x		x
376780		x	x								
376782		x	x								

371301 Blade Wrench

Handles and Accessories

REPLACEMENT COLLETS

HANDLES

371350 371394 371395 371410 372000 374310 374311 374320 374330 374340 375000 371305 371415

BLADES

	371350	371394	371395	371410	372000	374310	374311	374320	374330	374340	375000	371305	371415
376800	x					x	x	x	x	x		x	
376900	x					x	x	x	x	x		x	
376910	x					x	x	x	x	x		x	
376980		x	x										
376981		x	x										
376984		x	x										
377081		x	x										
377100	x					x	x	x	x	x		x	
377110	x					x	x	x	x	x		x	
377120	x					x	x	x	x	x		x	
377121	x					x	x	x	x	x		x	
377200	x					x	x	x	x	x		x	
377210	x					x	x	x	x	x		x	
377300	x					x	x	x	x	x		x	
377511	x					x	x	x	x	x			
377513	x					x	x	x	x	x			
377515	x					x	x	x	x	x			
377522	x					x	x	x	x	x			
377530	x					x	x	x	x	x			
377545	x					x	x	x	x	x			
377615													
377630													
378100				x									x
378300				x									x
378400				x									x
378600				x									x
378800				x									x
379071		x	x										
379081		x	x										
379082		x	x										
379100					x						x		
379400					x						x		
379500					x						x		
379700					x						x		
379710					x						x		

371301 Blade Wrench

Product Index

ITEM NO. DESCRIPTION CONTENTS PAGE

370110	Neuro-Sharp Blade, 15°, offset sickle	3	6
370120	Neuro-Sharp Blade, 15°, offset arachnoid	3	6
370410	Neuro-Sharp Blade, bayonet sickle	3	6
370420	Neuro-Sharp Blade, bayonet arachnoid	3	6
371301	Blade Wrench	1	14, 15
371305	Replacement Chuck Collet, stainless	1	13
371350	Handle, 13 cm L, .50 cm D	1	13
371360	Handle, 13 cm L, .66 cm D	1	13
371394	Arthro-Lok, 18 cm L, .40 cm handle	1	13, 14, 15
371395	Dynagrip, 21 cm L, .40 cm handle	1	13, 14, 15
371400	Bent Blade, 45°, sharp all around	4	7
371410	Handle, 10.2 cm L, .64 cm D	1	13
371415	Chuck Collet	1	13
371600	Bent Blade, 45°, sharp on one edge	4	7
371800	Bent Blade, 45°, sharp on two edges, blunt end	4	7
372000	BeaverTail, 17 cm L	1	13, 14, 15
374100	Ring Cutter	1	4
374150	Ring Cutter Saw Blade	1	4
374151	Ring Cutter Saw Blade, High Speed	1	4
374310	Handle, 10 cm L, .50 cm D	1	13
374311	Handle, 5 cm L, .50 cm D	1	13
374320	Handle, 15 cm L, .50 cm D	1	13
374330	Hex Handle, 10 cm L, .50 cm D	1	13
374340	Hex Handle, 15.5 cm L, .50 cm D	1	13
374410	Debakey Blade, 5 mm	24	8, 12
374451	Tympanoplasty, 2 mm W, angled 50°, bevel down	3	3

ITEM NO. DESCRIPTION CONTENTS PAGE

374562	Mini-Blade	24	7
374580	Arthro-Lok, 4 mm, angled 60°, hoe	6	9
374710	Debakey Blade, 8 mm	24	8, 12
374769	Mini-Knife, 180°, straight, sharp all around	3	2, 4, 10
374810	Debakey Blade, 10 mm	24	8, 12
374900	Blade, 270° cutting edge	6	12
375000	BeaverTail, 13 cm L	1	13, 14, 15
375910	Needle Blade	3	2, 6, 10
375920	Needle Blade	3	2, 6, 10
375930	Needle Blade	3	2, 6, 10
375940	Needle Blade	3	2, 6, 10
376083	Arthro-Lok, 4 mm, retrograde	6	9
376086	Arthro-Lok, 4 mm, pointed retrograde	6	9
376100	Mini-Blade	24	7
376300	Mini-Blade, pointed tip, sharp on two sides	6	2, 4, 10
376400	Mini-Blade, round tip, sharp on one side	24	2, 4, 10
376480	Arthro-Lok, 4 mm, round tip, sharp on one side, left	6	9
376481	Arthro Lok, 4 mm, round tip, sharp on one side, flat	6	9
376482	Arthro-Lok, 4 mm, round tip, sharp on one side, right	6	9
376500	Mini-Blade, smaller alternative to #11	24	2, 4, 10
376600	Lamellar Blade, 60°, bevel up	3	4, 10
376610	Lamellar Blade, sharp all around, smooth underside	3	4, 10
376700	Mini-Blade, smaller alternative to #15	24	2, 4, 10
376780	Arthro-Lok, 4 mm, pointed, curved L	6	9
376782	Arthro-Lok, 4 mm, pointed, curved R	6	9
376800	Tonsillectomy Blade, sickle shape	24	2


Product Index

ITEM NO. DESCRIPTION CONTENTS PAGE

376900	Mini-Blade, sharp all around	6	2, 4, 10
376910	Mini-Blade, 60°, cleft palate, beveled top & bottom	6	2, 10
376980	Arthro-Lok, 4 mm, mini-banana, curved	6	9
376981	Arthro-Lok, 4 mm, mini-banana, flat	6	9
376984	Arthro-Lok, 4 mm, banana	6	9
377081	Arthro-Lok, 4 mm, rosette	6	9
377100	Myringotomy Blade, 45°, lance tip, downward cut	6	3
377110	Myringotomy Blade, 45°, lance tip, upward cut	6	3
377120	Myringotomy Blade, 45°, spear tip	6	3
377121	Myringotomy Blade, 45°, lance tip, downward cut	6	3
377122	Myringotomy Knife, 45°, spear tip	3	3
377124	Myringotomy Knife, 45°, lance tip	3	3
377200	Tympanoplasty, 2.5 mm W, 60°, bevel down	3	3
377210	Tympanoplasty, 2.5 mm W, straight	3	3
377300	Sickle Edge, 84 mm L, straight	3	3, 6
377511	Micro-Sharp, 15°, 1.5 mm depth	6	5, 11
377512	Micro-Unitome, 15°, 1.5 mm depth	10	5, 11
377513	Micro-Sharp, 15°, 3 mm depth	6	5, 11
377514	Micro-Unitome, 15°, 3 mm depth	10	5, 11
377515	Micro-Sharp, 15°, 5 mm depth	6	5, 11
377516	Micro-Unitome, 15°, 5 mm depth	10	5, 11
377522	Micro-Sharp, 22.5°, 4 mm depth	6	5, 11
377523	Micro-Unitome, 22.5°, 4 mm depth	10	5, 11
377530	Micro-Sharp, 30°, 3.5 mm depth	6	5, 11
377531	Micro-Unitome, 30°, 3.5 mm depth	10	5, 11
377545	Micro-Sharp, 45°, 3 mm depth	6	5, 11

ITEM NO. DESCRIPTION CONTENTS PAGE

377546	Micro-Unitome, 45°, 3 mm depth	10	5, 11
377615	Micro-Blade, angled 15°	25	5, 11
377630	Micro-Blade, angled 30°	25	5, 11
378100	Chisel Blade, non-sterile	5	7
378300	Chisel Blade, non-sterile	5	7
378400	Chisel Blade, non-sterile	5	7
378600	Chisel Blade, non-sterile	5	7
378800	Chisel Blade, non-sterile	5	7
379071	Arthro-Lok, 4 mm, crescent, flat	6	9
379081	Arthro-Lok, 4 mm, mini-meniscus, flat	6	9
379082	Arthro-Lok, 4 mm, mini-meniscus, curved	6	9
379100	Cervical Biopsy Blade	4	12
379400	Chondroplastic Blade, sharp on two sides, curved	4	8
379500	Meniscus Blade, 7 mm, curved	4	8
379700	Meniscus Blade, 9 mm, curved	4	8
379710	Meniscus Blade, 9 mm, flat	4	8
581409	Bonn Forceps (0.12 tip)	10	12
581411	Tying Forceps	10	12
581417	Needle Holder	10	12
581422	Vannas Scissors	10	12

**BD Medical**

Ophthalmic Systems

1 Becton Drive
Franklin Lakes, NJ 07417
Tel: 800.237.2174 or 201.847.7142
Fax: 800.804.7489 or 201.847.4833
Email: special_services@BD.com

UK Customer Service

The Danby Building
Edmund Halley Road
Oxford Science Park
Oxford OX4 4DQ
England
Tel: 44.1.865.781666
Fax: 44.1.865.781627

Contact in other**European Countries**

11 rue Aristide Berges - BP4
38800 Le Pont de Claix
France
Tel: 33.4.76.68.32.27
Fax: 33.4.76.68.32.87