

Venous Products

BD Vacutainer® Accessories

BD Vacutainer® Accessories are designed for secure and safe specimen sampling. This product offering includes direct access sampling, transfer devices and tourniquets. These products offer convenience, ease of use, and complement the family of BD Vacutainer® and BD Microtainer® Products.



BD Vacutainer® Holders

BD Vacutainer® Holders and adapters provide quality, low-cost, single-use products that can help customers promote best practice infection control. The BD Vacutainer® One Use Holder is compatible with BD Vacutainer® Eclipse™ Blood Collection Needles, BD Vacutainer® Safety-Lok™ Blood Collection Sets, BD Vacutainer® Push Button Blood Collection Sets, and the BD Vacutainer® Multi Sample Luer Adapter.

Reference Number	Description	Packaging Box/Case Quantities
364815	BD Vacutainer® One Use Holder	250/1000
364597	BD Vacutainer® Ribbed Pediatric Tube Adapter	10/100
364888	BD Vacutainer® Holder	250/1000
368872	BD Pronto™ Quick Release Holder	20/100



BD Vacutainer® Blood Transfer Device

Transferring venous blood from a syringe and needle to a blood collection tube or blood culture bottle is a dangerous procedure. The BD Vacutainer® Blood Transfer Device was designed with healthcare workers' safety in mind. This single-use device reduces the risk of transfer-related injuries, while maintaining the specimen integrity required for accurate disease diagnosis. The BD Vacutainer® Blood Transfer Device offers safety without exception and maintains quality without compromise.



Reference Number	Description	Packaging Box/Case Quantities
364880	BD Vacutainer® Blood Transfer Device	200

BD Vacutainer® Luer-Lok™ Access Device

The BD Vacutainer® Luer-Lok™ Access Device is designed for sterile, secure, and safer specimen sampling. This device provides the security of a threaded, locking luer connection—the patented BD Vacutainer® Luer-Lok™, that replaces a luer slip device.

Reference Number	Description	Packaging Box/Case Quantities
364902	BD Vacutainer® Luer-Lok™ Access Device	200
366409	BD Vacutainer® Luer-Lok™ Access Device No Additive Plus Tube Kit: Luer-Lok™ Access Device, 6.0 mL, 13 x 100 Plus Plastic No Additive Tube with Clear BD Hemogard™ Closure	100
367300	BD Vacutainer® Luer Adapter	100/1000



BD Vacutainer® Stretch Latex-Free Tourniquet

Concerned about latex sensitivity? The BD Vacutainer® Stretch Latex-Free Tourniquet is ideal for practitioners and facilities that want to eliminate latex from their healthcare products. This latex-free tourniquet will not cause a latex-induced allergic reaction in latex-sensitive patients or employees. Convenient packaging allows for easy, one-at-a-time dispensing of tourniquets—encouraging a single-use policy that helps reduce the danger of cross-contamination between patients and healthcare workers.

Reference Number	Description	Packaging Box/Case Quantities
367204	BD Vacutainer® Stretch Latex-Free Tourniquet	25/500

