

BD Atomic Edge™ Accurate Depth Knives

Sharpness Redefined

Accurate Depth Knives for LRI and groove incisions

- Equal in performance to diamond blades
- Precise fixed blade exposure
- Brand new sterile knife with each case
- Single-use design



BD Atomic Edge™ Accurate Depth Knives

Cat. #	Description
370225	BD Atomic Edge™ Accurate Depth Knife, 250 micron
370230	BD Atomic Edge™ Accurate Depth Knife, 300 micron
370235	BD Atomic Edge™ Accurate Depth Knife, 350 micron
370240	BD Atomic Edge™ Accurate Depth Knife, 400 micron
370250	BD Atomic Edge™ Accurate Depth Knife, 500 micron
370255	BD Atomic Edge™ Accurate Depth Knife, 550 micron
370260	BD Atomic Edge™ Accurate Depth Knife, 600 micron
370262	BD Atomic Edge™ Accurate Depth Knife, 625 micron
370265	BD Atomic Edge™ Accurate Depth Knife, 650 micron

Also Available to Meet Your LRI Needs:

Cat. #	Description
372242	BD Beaver™ Swivel Fixation Ring
581409	BD Visitec™ Bonn Forceps (stainless steel)
581439	BD Visitec™ Bonn Forceps (plastic handle)

Additional information:

Incision Nomogram* developed by Robert M. Kershner, MD
(www.seeadoc.com)

Incision Nomogram* developed by James P. Gills, MD
(www.stlukeseye.com)

Incision Nomogram* developed by Louis D. "Skip" Nichamin, MD
(www.eyenz.com)

*These nomograms are provided to you as objective sources of clinically relevant information for your reference. None of the material contained in this information was generated by BD.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2009 BD 587906-0809



Helping all people
live healthy lives

For ordering information call:
1.866.906.8080 or 1.781.906.7900
or visit www.bd.com/safety

BD Medical
Ophthalmic Systems
411 Waverley Oaks Road
Waltham, MA 02452
Tel: 781.906.7900 or
866.906.8080
www.bd.com/ophthalmology