

Compatibility Chart

BD Pen Needles fit all diabetes pens and dosers in the U.S.



Drug Manufacturer	Pen Name	Drug	Drug NDC Number	Do BD Pen Needles Fit?*	Sample Bag of BD Pen Needles in Drug Starter Kits	Place Where Drug Manufacturers Reference BD Needles	Reference to BD Ultra-Fine™ Pen Needles
Sanofi-Aventis	Apidra® SoloStar®	Apidra®	0088-2502-05	Yes	N/A	Package Insert Patient Instructions for use Patient Toll Free Phone Number	"SoloStar® is compatible with all pen needles from Becton Dickinson and Company"
	Lantus® SoloStar®	Lantus®	0088-2220-60	Yes	N/A	Package Insert User Manual, Device Page, Website	"BD Ultra-Fine™ needles to be used in conjunction with SoloStar® and OptiClik® are sold separately and are manufactured by BD." "BD Ultra-Fine™ needles are compatible with SoloStar®. These are sold separately and are manufactured by BD."
	OptiClik® for Apidra®	Apidra® 3mL Cartridge	0088-2500-52	Yes	Yes	User Manual Patient Checklist	"Needles are available from BD Consumer Healthcare." "Prescription for BD pen needles."
	OptiClik® for Lantus®	Lantus® 3mL Cartridge	0088-2220-52	Yes	Yes	User Manual, Patient Checklist, Lantus® Package Insert	"Needles are available in different lengths and gauges from BD Consumer Healthcare."
Amylin Pharmaceuticals & Eli Lilly and Company	Byetta®	Byetta®	66780-0210-07 66780-0210-08	Yes	Yes		
Amylin Pharmaceuticals	SymlinPen™ 60 & 120	Symlin®	66780-155-02 66780-121-02	Yes	Yes, 5mm		
Eli Lilly and Company	HumaPen® Memoir™	Humalog® 3mL Cartridge	0002-7516-59	Yes	Yes, 8mm	User Manual, Website	"HumaPen® MEMOIR™ is suitable for use with Becton, Dickinson and Company insulin pen needles."
	HumaPen® LUXURA™ HD	Humalog® 3mL Cartridge	0002-7516-59	Yes	N/A, No Starter Kit Offered	User Manual, Website	"HumaPen® LUXURA™ HD is suitable for use with Becton, Dickinson and Company pen needles."
	Humalog® KwikPen™	Humalog® KwikPen™ Humalog® Mix 75/25® KwikPen™ Humalog® Mix 50/50 KwikPen™	0002-8799-59 0002-8798-59 0002-8797-59	Yes	Yes, 8mm	User Manual, Brochures, Prescribing Card	"When prescribing KwikPen™ prefilled with the Humalog brand of insulins you will need to write a separate prescription for BD needles."
	Humalog® Pen Humalog® Mix 50/50 Humalog® Mix 75/25® Pen	Humalog® Humalog® 50/50 Humalog® Mix 75/25®	0002-8725-59 0002-8739-59 0002-8794-59	Yes	Yes, 8mm	User Manual, Website	"This device is suitable for use with Becton, Dickinson and Company's insulin pen needles."
	Humulin® N Pen	Humulin® N	0002-8730-59	Yes	N/A, No Starter Kit Offered	User Manual, Website	"This device is suitable for use with Becton, Dickinson and Company's insulin pen needles."
	Humulin® 70/30 Pen	Humulin® 70/30	0002-8770-59	Yes	Humulin Starter Kit Offered	User Manual, Website	"This device is suitable for use with Becton, Dickinson and Company's insulin pen needles."
Novo Nordisk, Inc.	NovoLog® FlexPen®	NovoLog®	N/A	Yes			
	NovoPen® Junior	Novo Nordisk PenFill 3ml Insulin Cartridges	N/A	Yes			
	Novo Pen® 3	Novo Nordisk PenFill 3ml Insulin Cartridges	N/A	Yes			
	Levemir® FlexPen®	Levemir®	N/A	Yes			

* As of May, 2009

**Health Care Provider should select appropriate needle size/length

BD Ultra-Fine is a trademark of Becton, Dickinson and Company. All other brands are trademarks of their respective owners.

BD Pen Needles Fit All Pens and dosers sold in the U. S.* including the new Apidra® SoloStar®!



Apidra® SoloStar®



Lantus® SoloStar®



OptiClik® for Lantus® and Apidra®



SymlinPen® 120 (pramlintide acetate) pen injector



SymlinPen® 60 (pramlintide acetate) pen injector



Byetta® exenatide injection 10mcg



Byetta® exenatide injection 5mcg



Levemir® FlexPen®



Novolog® FlexPen®
Novolog® Mix 70/30 FlexPen®



NovoPen® Junior



NovoPen® 3



HumaLog® KwikPen™
HumaLog® Mix 75/25® KwikPen™
HumaLog® Mix 50/50 KwikPen™



Humalog® Pen
Humalog® Mix 75/25® Pen
Humalog® Mix 50/50 Pen



Humulin® N Pen
Humulin® 70/30 Pen



HumaPen® LUXURA™ HD



HumaPen® MEMOIR™



AutoPen®



Forteo® teriparatide (rDNA origin) injection



Novolin®
InnoLet®

Apidra® SoloStar® is compatible with all
BD Ultra-Fine™ Pen Needles



Apidra® SoloStar®



BD

Helping all people
live healthy lives

The BD Ultra-Fine™ Mini Pen Needle

5 mm (3/16") x 31G


NDC number: 08290-3201-19



The BD Ultra-Fine™ Short Pen Needle

8 mm (5/16") x 31G

NDC number: 08290-3201-09



The BD Ultra-Fine™ Original Pen Needle

12.7 mm (1/2") x 29G

NDC number: 08290-3282-03



For more information call 1-888-232-2737
or visit www.BDdiabetes.com/us

* As of May, 2009

BD, BD Logo and BD Ultra-Fine are trademarks of Becton, Dickinson and Company. 2009 BD.
All other brands are trademarks of their respective owners.
09D06512141STAR

Needle is not actual size